

Medtronic

Getting people back to fuller lives

Peripheral Vascular Health
U.S. product catalog



Learn more and
connect with us at
medtronic.com/PVH





Peripheral Vascular Health
is dedicated to the treatment
of peripheral arterial and
venous diseases.

Our goal is to help physicians save
limbs, reduce pain, and treat disease
– so patients can get back to
enjoying what they love.

Table of contents



Directional Atherectomy Systems

HawkOne™ Directional Atherectomy System.....	4
SilverHawk™ Peripheral Plaque Excision System	4

Drug-Coated Balloons

IN.PACT™ Admiral™ Drug-Coated Balloon	6
IN.PACT™ 018 Drug-Coated Balloon	7
IN.PACT™ AV Drug-Coated Balloon	8

Endovascular Arteriovenous Fistula (EndoAVF)

Ellipsys™ Vascular Access System	10
--	----

Embolic Protection

SpiderFX™ Embolic Protection Device	12
Mo.Ma Ultra™ Proximal Cerebral Protection Device	12

CTO Devices

Viance™ Crossing Catheter.....	14
Enteer™ Re-entry System.....	14

Embolization

Concerto™ Detachable Coil System.....	16
MVP™ Micro Vascular Plug System.....	17
Rebar™ Reinforced Micro Catheter	17

Superficial Vein Products

VenaSeal™ Closure System.....	19
ClosureFast™ Radiofrequency Ablation System	20

Stents

Self-expanding Stent Systems

EverFlex™ Self-expanding Peripheral Stent with Entrust™ Delivery System.....	22
EverFlex™ Self-expanding Peripheral Stent System	23
Protégé™ EverFlex™ Self-expanding Biliary Stent System	24
Protégé™ GPS™ Self-expanding Peripheral and Biliary Stent System	25
Protégé™ RX Self-expanding Carotid Stent System	26

Balloon-expandable Stent System

Visi-Pro™ Balloon-expandable Peripheral Stent System	27
--	----

Table of contents



Unmounted Stents

IntraStent™ DoubleStrut LD Biliary Stent	28
IntraStent™ Mega LD Biliary Stent	28
IntraStent™ Max LD Biliary Stent	28

Venous Stents

Abre™ Venous Self-expanding Stent System	29
--	----

PTA Balloons

EverCross™ 0.035" OTW PTA Dilatation Catheter	31
Admiral Xtreme™ PTA Balloon Catheter	32
Fortrex™ 0.035" OTW HP PTA Balloon Catheter	33
Pacific™ Plus PTA Catheter OTW 0.018"	34
NanoCross™ Elite 0.014" OTW PTA Balloon Catheter	36
RapidCross™ PTA Rapid Exchange Balloon Dilatation Catheter	38

PTA Specialty Balloon

Chocolate™** PTA Balloon	40
--------------------------------	----

Infusion Therapy

Cragg-McNamara™ Valved Infusion Catheters	42
---	----

Snares

Amplatz Goose Neck™ Snare Kit	44
Amplatz Goose Neck™ MicroSnare Kit	44

Guidewires

Wholey™** Guidewire System	46
Babywire™ Double-ended Nitinol Guidewire	46
Nitrex™ Guidewire	47


Support Catheters and Guide Catheters

TrailBlazer™ Support Catheter	49
TrailBlazer™ Angled Support Catheter	49
Launcher™ Peripheral Guide Catheter	50

Brief statements and reference statements	51
---	----



Directional Atherectomy Systems



HawkOne™
Directional Atherectomy
System

HawkOne™

Directional Atherectomy System

The versatile HawkOne device – the latest addition to the Medtronic directional atherectomy portfolio – restores blood flow in patients with peripheral arterial disease (PAD) by removing plaque from calcified arteries. The HawkOne device is a comprehensive system that treats all plaque morphologies, from soft to severely calcified, above and below the knee, and offers physicians the flexibility to create a channel or maximize luminal gain.

Model name	Reference number	Vessel diameter (mm)	Sheath compatibility (F)	Crossing profile (mm)	Working length [†] (cm)	Effective length [‡] (cm)	Tip length (cm)	Max. cut length (mm)	Packing device
HawkOne LS Standard Tip	H1-LS	3.5 to 7.0	7	2.6	114	107	6.6	50	■
HawkOne LX Extended Tip	H1-LX	3.5 to 7.0	7	2.6	114	104	9.6	75	■
HawkOne M Standard Tip	H1-M	3.0 to 7.0	6	2.2	135	129	5.9	40	■
HawkOne S Standard Tip	H1-S	2.0 to 4.0	6	2.2	151	145	5.9	40	■

Max guidewire is 0.014" for all HawkOne devices.

[†]Working length – distal end of preloaded flush tool, in the proximal position, to the distal end of tip.

[‡]Effective length – distal end of preloaded flush tool, in the proximal position, to the proximal end of cutter window.

SilverHawk™

Peripheral Plaque Excision System

The first-generation SilverHawk plaque excision system treats peripheral arterial disease (PAD) by removing soft-to-mild plaque buildup in arteries below the knee. The SilverHawk system technology uses a directional cutting blade to shave plaque from the vessel. The plaque is captured in the nosecone and safely removed from the vessel.


Product name	Reference number	Vessel diameter (mm)	Sheath compatibility (F)	Crossing profile (mm)	Working length [†] (cm)	Effective length [‡] (cm)	Tip length (cm)	Max cut length (mm)	Packing device
EXL	P4044	2.0 to 3.0	6	2.0	135	129	6.0	15	■
DS	P4028	1.5 to 2.0	6	1.9	135	132	2.6	10	

Max guidewire is 0.014" for all SilverHawk devices.

[†]Working length – distal end of strain relief to the distal end of tip.

[‡]Effective length – distal end of strain relief to the proximal end of cutter window.

Drug-Coated Balloons



IN.PACT™ AV
Drug-Coated Balloon



IN.PACT™ Admiral™
Drug-Coated Balloon



IN.PACT™ 018
Drug-Coated Balloon

IN.PACT™ Admiral™

Drug-Coated Balloon



The IN.PACT Admiral drug-coated balloon (DCB) is a clinically proven, primary endovascular therapy that treats femoropopliteal disease, reduces interventions, and preserves future treatment options.¹

Technical specifications

Paclitaxel drug dose	3.5 µg/mm ²	Catheter design	Over the wire (OTW)
Excipient	Urea	Catheter lengths	80, 130 cm
Balloon diameters	4.0–7.0 mm	Guidewire compatibility	0.035"
Balloon lengths	40, 60, 80, 120, 150, 200, 250 mm [†]	Nominal balloon pressure	8 atm: 40, 60, 80, 120, and 150 mm 5 atm: 200 and 250 mm
Balloon fold configuration	4.0 mm: 3 folds 5.0, 6.0, and 7.0 mm: 6 folds		

[†]120, 150, 200, and 250 mm lengths are not offered on the 7.0 mm diameter balloon.

¹ Laird JA, Schneider PA, Jaff MR, et al. Long-Term Clinical Effectiveness of a Drug-Coated Balloon for the Treatment of Femoropopliteal Lesions. *Circ Cardiovasc Interv.* June 2019;12(6):e007702.

Reference number usable length 80 cm	Reference number usable length 130 cm	Balloon diameter (mm)	Balloon length (mm)	Recommended introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
ADM04004008P	ADM04004013P	4.0	40	5	8	14
ADM04006008P	ADM04006013P	4.0	60	5	8	14
ADM04008008P	ADM04008013P	4.0	80	5	8	14
ADM04012008P	ADM04012013P	4.0	120	5	8	14
ADM04015008P	ADM04015013P	4.0	150	5	8	14
-	ADM04020013P	4.0	200	5	5	11
-	ADM04025013P	4.0	250	5	5	11
ADM05004008P	ADM05004013P	5.0	40	6	8	14
ADM05006008P	ADM05006013P	5.0	60	6	8	14
ADM05008008P	ADM05008013P	5.0	80	6	8	14
ADM05012008P	ADM05012013P	5.0	120	6	8	14
ADM05015008P	ADM05015013P	5.0	150	6	8	14
-	ADM05020013P	5.0	200	6	5	11
-	ADM05025013P	5.0	250	6	5	11
ADM06004008P	ADM06004013P	6.0	40	6	8	14
ADM06006008P	ADM06006013P	6.0	60	6	8	14
ADM06008008P	ADM06008013P	6.0	80	6	8	14
ADM06012008P	ADM06012013P	6.0	120	6	8	14
ADM06015008P	ADM06015013P	6.0	150	6	8	14
-	ADM06020013P	6.0	200	6	5	11
-	ADM06025013P	6.0	250	6	5	11
ADM07004008P	ADM07004013P	7.0	40	7	8	14
ADM07006008P	ADM07006013P	7.0	60	7	8	14
ADM07008008P	ADM07008013P	7.0	80	7	8	14

Risks of the IN.PACT Admiral DCB may include access site pain, hemorrhage, local infection at access site, local or distal embolic events, perforation or rupture of the artery, amputation/loss of limb, and death.

IN.PACT™ 018

Drug-Coated Balloon



The IN.PACT 018 drug-coated balloon (DCB) is a primary endovascular therapy for the treatment of femoropopliteal disease. It features a low-profile design engineered to cross tight lesions and allows the option to treat via femoral or radial access with 130 cm and 200 cm catheter lengths.

Technical specifications

Paclitaxel drug dose	3.5 µg/mm ²	Catheter design	Over the wire (OTW)
Excipient	Urea	Catheter lengths	130, 200 cm
Balloon diameters	4.0–7.0 mm	Guidewire compatibility	0.018"
Balloon lengths	40, 60, 80, 100, 120, 150 mm [†]	Nominal balloon pressure	8 atm
Balloon fold configuration	4.0 mm: 3 folds 5.0, 6.0, and 7.0 mm: 6 folds		

[†]100, 120 and 150 mm lengths are not offered on the 7.0 mm diameter balloon.

Ref. number usable length 130 cm	Ref. number usable length 200 cm	Balloon diameter (mm)	Balloon length (mm)	Recommended introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
IPU04004013P	IPU04004020P	4.0	40	5	8	10
IPU04006013P	IPU04006020P	4.0	60	5	8	10
IPU04008013P	IPU04008020P	4.0	80	5	8	10
IPU04010013P	IPU04010020P	4.0	100	5	8	10
IPU04012013P	IPU04012020P	4.0	120	5	8	10
IPU04015013P	IPU04015020P	4.0	150	5	8	10
IPU05004013P	IPU05004020P	5.0	40	5	8	10
IPU05006013P	IPU05006020P	5.0	60	5	8	10
IPU05008013P	IPU05008020P	5.0	80	5	8	10
IPU05010013P	IPU05010020P	5.0	100	5	8	10
IPU05012013P	IPU05012020P	5.0	120	5	8	10
IPU05015013P	IPU05015020P	5.0	150	5	8	10
IPU06004013P	IPU06004020P	6.0	40	5	8	10
IPU06006013P	IPU06006020P	6.0	60	5	8	10
IPU06008013P	IPU06008020P	6.0	80	5	8	10
IPU06010013P	IPU06010020P	6.0	100	5	8	10
IPU06012013P	IPU06012020P	6.0	120	5	8	10
IPU06015013P	IPU06015020P	6.0	150	5	8	10
IPU07004013P	IPU07004020P	7.0	40	6	8	10
IPU07006013P	IPU07006020P	7.0	60	6	8	10
IPU07008013P	IPU07008020P	7.0	80	6	8	10

Risks of IN.PACT 018 DCB may include, but are not limited to, abrupt vessel closure, vessel spasm, perforation or rupture of the artery, dissection, pseudoaneurysm, hematoma, thrombosis, and stroke.

IN.PACT™ AV

Drug-Coated Balloon



The IN.PACT AV drug-coated balloon (DCB) is a clinically demonstrated endovascular therapy for AV fistula maintenance for patients with end-stage renal disease (ESRD). It delivers an antiproliferative drug (paclitaxel) to the vessel to inhibit AV fistula stenosis. The IN.PACT AV DCB may enable dramatically fewer AV fistula reinterventions, which could keep patients out of the hospital longer.¹ It can make an impact clinically, financially,² and emotionally.

¹ Lookstein RA, Haruguchi H, Ouriel K, et al. Drug-Coated Balloons for Dysfunctional Dialysis Arteriovenous Fistulas. *N Engl J Med*. August 20, 2020;383(8):733-742.

² Thamer M, Lee TC, Wasse H, et al. Medicare Costs Associated With Arteriovenous Fistulas Among US Hemodialysis Patients. *Am J Kidney Dis*. July 2018;72(1):10-18.

Ref. number usable length 40 cm	Ref. number usable length 80 cm	Ref. number usable length 130 cm	Balloon diameter (mm)	Balloon length (mm)	Recommended introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
IAV04004004P	IAV04004008P	-	4.0	40	5	8	14
IAV04006004P	IAV04006008P	-	4.0	60	5	8	14
IAV04008004P	IAV04008008P	-	4.0	80	5	8	14
IAV04012004P	IAV04012008P	-	4.0	120	5	8	14
IAV05004004P	IAV05004008P	-	5.0	40	6	8	14
IAV05006004P	IAV05006008P	-	5.0	60	6	8	14
IAV05008004P	IAV05008008P	-	5.0	80	6	8	14
IAV05012004P	IAV05012008P	-	5.0	120	6	8	14
IAV06004004P	IAV06004008P	-	6.0	40	6	8	14
IAV06006004P	IAV06006008P	-	6.0	60	6	8	14
IAV06008004P	IAV06008008P	-	6.0	80	6	8	14
IAV06012004P	IAV06012008P	-	6.0	120	6	8	14
IAV07004004P	IAV07004008P	-	7.0	40	7	8	14
IAV07006004P	IAV07006008P	-	7.0	60	7	8	14
IAV07008004P	IAV07008008P	-	7.0	80	7	8	14
IAV08004004P	IAV08004008P	IAV08004013P	8.0	40	7	8	10
IAV08006004P	IAV08006008P	IAV08006013P	8.0	60	7	8	10
IAV08008004P	IAV08008008P	IAV08008013P	8.0	80	7	8	10
IAV09004004P	IAV09004008P	IAV09004013P	9.0	40	7	8	10
IAV09006004P	IAV09006008P	IAV09006013P	9.0	60	7	8	10
IAV09008004P	IAV09008008P	IAV09008013P	9.0	80	7	8	10
IAV10004004P	IAV10004008P	IAV10004013P	10.0	40	7	6	9
IAV12004004P	IAV12004008P	IAV12004013P	12.0	40	9	6	9

EndoAVF



Ellipsys™
Vascular Access System

Ellipsys™

Vascular Access System



Designed for end-stage renal disease (ESRD) patients requiring hemodialysis, the Ellipsys system is a unique single-catheter, nonsurgical option for physicians to create an arteriovenous (AV) fistula. It transforms the standard surgical AV fistula creation into a reproducible, minimally invasive procedure that requires no implant or suture, and allows patients to leave with just a single needle stick.¹

¹ Hull JE, Jennings WC, Cooper RI, et al. The Pivotal Multicenter Trial of Ultrasound-Guided Percutaneous Arteriovenous Fistula Creation for Hemodialysis Access. J Vasc Interv Radiol. February 2018;29(2):149-158.e5.

Reference number	Description	Recommended introducer sheath (F)	Guidewire compatibility (in)	Vessel size (mm)
AMI6005	Ellipsys catheter†	6‡	0.014	2.0 +

†Requires AMI-1001 Ellipsys™ power controller (110-240v, 50/60HZ, reusable).

‡Maximum sheath length 10 cm. Prior to using a specific sheath, ask your sales or clinical representative if the sheath is compatible with the Ellipsys catheter. Sheaths that are too long may cover the working portion at the device tip.

Ellipsys system

Reference number	Description	Units per pack
AMI1001	Ellipsys power controller (reusable)	1 each
K21-00022	EndoAVF access kit hydrophilic sheath§	5 each

§Includes 21 ga x 7 cm echo-tip micropuncture needle, 0.018" guidewire, and 6F x 10 cm thin-walled radial artery sheath.

Ellipsys vascular access system accessories


Reference number	Description	Units per pack
N140801	Nitrex™ guidewire, diameter 0.014", length 80 cm	3 each

Risks of the Ellipsys system may include total/partial occlusion or stenosis of the anastomosis, failure to achieve fistula maturation, Steal Syndrome, hematoma, infection, and need for vessel superficialization or other maturation assistance procedures.

Embololic Protection



Mo.Ma Ultra™
Proximal Cerebral Protection Device



SpiderFX™
Embololic Protection Device

SpiderFX™

Embolic Protection Device

The SpiderFX device is used to capture and remove debris that becomes dislodged during an interventional procedure. The SpiderFX device can be delivered over any 0.014" or 0.018" guidewire, or through any 0.035" guidewire-compatible catheter.† It is indicated for use in carotid arteries, coronary saphenous vein bypass grafts, and lower extremity procedures.

†Lower extremity procedures.

Reference number	Filter size (mm)	Target vessel size (mm)	Capture wire		Delivery end	Recovery end	Guide catheter/sheath
			Wire length OTW/RX (cm)	Wire diameter (in/mm)	Crossing profile (F)	Diameter (F)	Minimum ID (in)
SPD2-US-030-320	3.0	3.0 SVG & Carotid	320/190	0.014/0.36	3.2	4.2	0.066
SPD2-US-040-320	4.0	3.1–4.0 SVG & Carotid	320/190	0.014/0.36	3.2	4.2	0.066
SPD2-US-050-320	5.0	4.1–5.0 SVG & Carotid 3.0–4.0 Lower Extremity	320/190	0.014/0.36	3.2	4.2	0.066
SPD2-US-060-320	6.0	4.5–6.0 SVG & Carotid 3.5–5.0 Lower Extremity	320/190	0.014/0.36	3.2	4.2	0.066
SPD2-US-070-320	7.0	5.5–6.0 SVG 5.5–7.0 Carotid 4.5–6.0 Lower Extremity	320/190	0.014/0.36	3.2	4.2	0.066

The 320 cm wire lengths are scored to allow for snapping down to a 190 cm wire length if desired.

Mo.Ma Ultra™

Proximal Cerebral Protection Device


The Mo.Ma Ultra proximal cerebral protection device is used to contain and remove debris that can dislodge during interventional procedures in the carotid arteries. The Mo.Ma Ultra device with double-occlusion balloon system allows for proximal embolic protection to be established prior to crossing a carotid lesion.

Technical specifications

Balloon material	Compliant elastomeric rubber	
Balloon marker distance	6 cm	
Recommended guidewire	0.035" (0.89 mm)	
Balloon occlusion range	5–13 mm diameter (CCA prox. balloon)	3–6 mm diameter (ECA dist. balloon)

Reference number	Minimum sheath size (F)	Inner diameter of the working channel
MUS0130069X6	9	0.083"

CTO Devices



Enteer™
Re-entry System



Viance™
Crossing Catheter

Viance™

Crossing Catheter

The Viance crossing catheter is designed to efficiently cross chronic total occlusions via the true lumen. The low-profile catheter, with its fast-spin torque handle, is designed to find small microchannels in a lesion, while leaving the control of crossing in the physician's hands.

Reference number	Description	Working length (cm)	Guidewire compatibility (in)	Crossing profile (max in)	Sheath compatibility
VNC-FX-150	Flexible	150	0.014	0.038	5 F
VNC-SD-150	Standard	150	0.014	0.038	5 F

Enteer™ Re-entry System

The Enteer system, consisting of the Enteer re-entry balloon catheter and the Enteer guidewire, provides the physician with control to reliably target the true lumen from the subintimal channel above or below the knee. The Enteer catheter's unique balloon design self-oriens toward the true lumen within the subintimal space when inflated. Offset exit ports are located on either side of the device, allowing the Enteer guidewire to reenter into the true lumen.


Re-entry Balloon Catheter

Reference number		Balloon size (W x H x L mm)	Working length (cm)	Guidewire compatibility (in)	Crossing profile (max in)	Sheath compatibility
ENB-375-20-135	ATK	3.75 x 1.5 x 20	135	≤ 0.018	0.066	5 F
ENB-275-20-150	BTK	2.75 x 1.0 x 20	150	≤ 0.018	0.066	5 F


Guidewire

Reference number	Description	Diameter (in)	Length (cm)	Tip reach (mm)
ENW-FX-014-300	Flexible	0.014	300	1.5
ENW-SD-014-300	Standard	0.014	300	1.5
ENW-SF-014-300	Stiff	0.014	300	2.5

Embolization



Rebar™
Reinforced
Micro Catheter



MVP™
Micro Vascular
Plug System



Concerto™
3D Detachable
Coil System

Concerto™
Helix Detachable
Coil System

Concerto™

Detachable Coil System



Treat more patients, treat them safely, and treat them faster with Concerto coils. The combination of complex and helical coils enables you to treat multiple anatomies while the coils' softness helps smoothly reach distal and challenging vessels.¹ The detachment system has a 99.7% success rate,¹ and fibers reduce occlusion time, enhancing the mechanical action of bare coils.²

¹Test data on file at Medtronic in reports TR08-134, and TR07-084. Bench test results may not be indicative of clinical performance.

²Girdhar G, Read M, Sohn J, et al. In-vitro thrombogenicity assessment of polymer filament modified and native platinum embolic coils. *J Neurol Sci.* April 15, 2014;339(1-2):97-101.

Order number	Description	Diameter (mm)	Length (cm)	Shape	Min. micro catheter compatibility (in)
NV-2-4-Helix	Concerto nylon helical	2	4	Helix	0.0165
NV-2-6-Helix	Concerto nylon helical	2	6	Helix	0.0165
NV-2-8-Helix	Concerto nylon helical	2	8	Helix	0.0165
NV-3-4-Helix	Concerto nylon helical	3	4	Helix	0.0165
NV-3-8-Helix	Concerto nylon helical	3	8	Helix	0.0165
NV-4-8-Helix	Concerto nylon helical	4	8	Helix	0.0165
NV-4-10-Helix	Concerto nylon helical	4	10	Helix	0.0165
NV-5-15-Helix	Concerto nylon helical	5	15	Helix	0.021
NV-5-20-Helix	Concerto nylon helical	5	20	Helix	0.021
NV-6-20-Helix	Concerto nylon helical	6	20	Helix	0.021
NV-7-30-Helix	Concerto nylon helical	7	30	Helix	0.021
NV-8-30-Helix	Concerto nylon helical	8	30	Helix	0.021
NV-9-30-Helix	Concerto nylon helical	9	30	Helix	0.021
NV-10-30-Helix	Concerto nylon helical	10	30	Helix	0.021
PV-12-30-Helix	Concerto PGLA helical	12	30	Helix	0.021
PV-14-30-Helix	Concerto PGLA helical	14	30	Helix	0.021
PV-16-40-Helix	Concerto PGLA helical	16	40	Helix	0.021
PV-18-40-Helix	Concerto PGLA helical	18	40	Helix	0.021
PV-20-50-Helix	Concerto PGLA helical	20	50	Helix	0.021
PV-2-2-3D	Concerto PGLA 3D	2	2	3D	0.0165
PV-2-4-3D	Concerto PGLA 3D	2	4	3D	0.0165
PV-2-6-3D	Concerto PGLA 3D	2	6	3D	0.0165
PV-3-4-3D	Concerto PGLA 3D	3	4	3D	0.0165
PV-3-6-3D	Concerto PGLA 3D	3	6	3D	0.0165
PV-3-8-3D	Concerto PGLA 3D	3	8	3D	0.0165
PV-4-8-3D	Concerto PGLA 3D	4	8	3D	0.0165
PV-4-10-3D	Concerto PGLA 3D	4	10	3D	0.0165
PV-4-12-3D	Concerto PGLA 3D	4	12	3D	0.0165
PV-5-15-3D	Concerto PGLA 3D	5	15	3D	0.0165
PV-6-20-3D	Concerto PGLA 3D	6	20	3D	0.0165
PV-7-30-3D	Concerto PGLA 3D	7	30	3D	0.0165
PV-8-30-3D	Concerto PGLA 3D	8	30	3D	0.0165
PV-9-30-3D	Concerto PGLA 3D	9	30	3D	0.0165
PV-10-30-3D	Concerto PGLA 3D	10	30	3D	0.0165
PV-12-40-3D	Concerto PGLA 3D	12	40	3D	0.021
PV-14-40-3D	Concerto PGLA 3D	14	40	3D	0.021
PV-16-40-3D	Concerto PGLA 3D	16	40	3D	0.021
PV-18-40-3D	Concerto PGLA 3D	18	40	3D	0.021

Risks of the Concerto system may include puncture site hematoma, thromboembolic episodes, vessel perforation, neurological deficits, vasospasms, vascular thrombosis, hemorrhage, and ischemia.

MVP™

Micro Vascular Plug System

The MVP micro vascular plug system is designed to occlude vessels in the peripheral vasculature quickly and predictably. Micro catheter deliverability enables super-selective treatment of distal vessels. Further, the device is fully resheathable to facilitate accurate placement.

Order number	Description	Recommended vessel size (mm)	Outer diameter unconstrained (mm)	Length unconstrained (mm)	Delivery wire length (cm)	Min. recommended microcatheter ID (in)	Max. recommended length of microcatheter/delivery catheter (cm)
MVP-3Q	MVP-3Q	1.5-3	5.3	12	180	0.021	153
MVP-5Q	MVP-5Q	3-5	6.5	12	180	0.027	153
MVP-7Q	MVP-7Q	5-7	9.2	16	165	0.041 (4 Fr)	120
MVP-9Q	MVP-9Q	7-9	13	18	165	0.043 (5 Fr)	120

Rebar™

Reinforced Micro Catheter

The Rebar micro catheter is an end-hole, single-lumen, stainless steel reinforced catheter designed to be introduced over a steerable guidewire into the vasculature. It features a lubricious hydrophilic outer coating and steam-shapeable tip.

Order number	Product	Proximal OD/distal OD (F)	Distal ID (in)	Total length (cm)	Usable length (cm)	Maximum guidewire (in)
105-5081-130	Rebar-18	2.7/2.4	0.021	137	130	0.018
105-5083-153	Rebar-18	2.7/2.4	0.021	160	153	0.018
105-5082-145	Rebar-27	2.8/2.8	0.027	150	145	0.021

Superficial Venous Products



VenaSeal™
Closure System



ClosureFast™
Radiofrequency Ablation System

VenaSeal™

Closure System



The VenaSeal closure system offers relief for patients suffering from venous reflux disease by using a specially formulated medical adhesive to permanently close the vein. With no need for heat, the VenaSeal procedure delivers a comfortable patient experience^{1,2} and immediate and lasting vein closure with a demonstrated 94.6% closure rate at five years.^{1,2}

¹ Proebstle TM. The European Multicenter Study on Cyanoacrylate Embolization of Refluxing Great Saphenous Veins without Tumescence Anesthesia and without Compression Therapy. Results presented at Charing Cross 2016; London, UK.

² Almeida JI, Javier JJ, Mackay EG, et al. Thirty-sixth-month follow-up of first-in-human use of cyanoacrylate adhesive for treatment of saphenous vein incompetence. *J Vasc Surg Venous Lymphat Disord*. September 2017;5(5):658-666.

Description	Specs
VenaSeal Adhesive	5 cc of the VenaSeal adhesive (a specially formulated n-butyl-2-cyanoacrylate) is contained within a screw-capped vial.
Syringe	The 3 cc syringes are graduated Monoject™ luer lock syringes, each with a standard threaded luer lock connector.
Introducer	The introducer is 7 F with an effective length of 80 cm, with 10 mm spaced circumferential markings along its length for measuring retraction length during the VenaSeal procedure.
Guidewire	The guidewire is a 0.035 in, 180 cm straight floppy-tip guidewire.
Dispenser Tips	The dispenser tips are each comprised of a stainless steel, 1.5 mm inner diameter, 3.8 cm length hypotube with a luer lock connector.
Dispenser Gun	The dispenser gun consists of a pistol-type ergonomic handle with an integrated barrel and trigger. Each depression of the trigger delivers a controlled 0.10 cc (range: 0.06–0.12 cc) amount of adhesive.
Dilator	The dilator is 5 F with an effective length of 87 cm.
Catheter	The catheter is 5 F with an effective length of 91 cm, laser markings at 3 cm and 85 cm from the tip, and high echogenic visibility.

ClosureFast™

Radiofrequency Ablation (RFA) System

Deliver lasting results with gap-free thermal ablation using the ClosureFast RFA system. The segmental procedure offers proven results with positive patient outcomes and experience.¹⁻³ It is the only radiofrequency ablation system with published long-term clinical data demonstrating safety and efficacy, with a 91.9% closure rate at five years.²

¹ Almeida JJ, Kaufman J, Göckeritz O, et al. Radiofrequency endovenous ClosureFAST versus laser ablation for the treatment of great saphenous reflux: a multicenter, single-blinded, randomized study (RECOVERY study). *J Vasc Interv Radiol.* June 2009;20(6):752-759.

² Proebstle TM, Alm BJ, Göckeritz O, et al. Five-year results from the prospective European multicentre cohort study on radiofrequency segmental thermal ablation for incompetent great saphenous veins. *Br J Surg.* February 2015;102(3): 212-218.

³ Rasmussen LH, Lawaetz M, Bjoern L, et al. Randomized clinical trial comparing endovenous laser ablation, radiofrequency ablation, foam sclerotherapy and surgical stripping for great saphenous varicose veins. *Br J Surg.* August 2011;98(8):1079-1087.

Catheter model	CF7-7-60	CF7-7-100	CF7-3-60
Compatible Guidewire	0.025 in	0.025 in	0.025 in
Default Target Temperature Setting	120°C	120°C	120°C
Heating Element Diameter	7 F (2.3 mm)	7 F (2.3 mm)	7 F (2.3 mm)
Heating Element Length	7 cm	7 cm	3 cm
Heating Element Overlap	0.5 cm	0.5 cm	0.5 cm
Insertable Length	60 cm	100 cm	60 cm
Introducer Sheath (minimum ID size)	7 F (2.3 mm)	7 F (2.3 mm)	7 F (2.3 mm)
Maximum Power Setting	40 W	40 W	18 W
Software Version: RFG2	4.0.0 or higher	4.0.0 or higher	4.0.0 or higher
Software Version: RFG3	1.11.0 or higher	1.11.0 or higher	1.11.0 or higher

ClosureRFS™

Endovenous Radiofrequency Ablation (RFA) Stylet

The ClosureRFS stylet RFA device is indicated for the treatment of incompetent perforator veins for patients suffering from venous reflux disease.

Catheter model	Product names	French size	Working length	Compatible guidewire	Max IV catheter
RFS2-6-12	ClosureRFS stylet	6 F (2.0 mm)	12 cm	0.035 in	12 gauge

ClosureRFG™


Radiofrequency Generator

The ClosureRFG radiofrequency generator supplies and controls the RF energy delivered to the ClosureFast catheter and ClosureRFS stylet, providing consistent and controlled treatment of venous reflux disease.


Product catalog number	Voltage	Height	Width	Depth	Weight
RFG3	Universal (100-240V)	10.5 in (26.7 cm)	13.4 in (34 cm)	6.8 in (17.3 cm)	15 lb (6.8 kg) max

Risks of the ClosureFast procedure may include hematoma, phlebitis, skin injury, nerve injury, thrombophlebitis, thrombosis, and/or pulmonary embolism.

Stents



EverFlex™
Self-expanding Peripheral
Stent System



Abre™
Venous Self-expanding
Stent System

Visi-Pro™
Balloon-expandable Peripheral
Stent System



EverFlex™ Self-expanding Peripheral Stent with Entrust™ Delivery System

You asked for simple deployment with reduced variability – and we delivered. When you need to stent, trust the precision, strength, and flexibility of the #1 peripheral stent†: EverFlex self-expanding peripheral stent with Entrust delivery system. The Entrust system is a one-handed stent delivery system with a low 5 F profile. The device is indicated for use in the superficial femoral artery and/or proximal popliteal arteries, and is available in lengths from 20 mm to 150 mm.

†EverFlex Self-expanding Peripheral Stent, U.S. only. DRG market share data for peripheral self-expanding bare metal stents.

Catheter length			Stent size		Compatibility		
80 cm reference number	120 cm reference number	150 cm reference number	Diameter (mm)	Length (mm)	Recomm. introducer sheath (F)	Guidewire (in)	Recomm. lumen size (mm)
EVD35-06-020-080	EVD35-06-020-120	EVD35-06-020-150	6	20	5	0.035	4.5-5.5
EVD35-06-040-080	EVD35-06-040-120	EVD35-06-040-150	6	40	5	0.035	4.5-5.5
EVD35-06-060-080	EVD35-06-060-120	EVD35-06-060-150	6	60	5	0.035	4.5-5.5
EVD35-06-080-080	EVD35-06-080-120	EVD35-06-080-150	6	80	5	0.035	4.5-5.5
EVD35-06-100-080	EVD35-06-100-120	EVD35-06-100-150	6	100	5	0.035	4.5-5.5
EVD35-06-120-080	EVD35-06-120-120	EVD35-06-120-150	6	120	5	0.035	4.5-5.5
EVD35-06-150-080	EVD35-06-150-120	EVD35-06-150-150	6	150	5	0.035	4.5-5.5
EVD35-07-020-080	EVD35-07-020-120	EVD35-07-020-150	7	20	5	0.035	5.5-6.5
EVD35-07-040-080	EVD35-07-040-120	EVD35-07-040-150	7	40	5	0.035	5.5-6.5
EVD35-07-060-080	EVD35-07-060-120	EVD35-07-060-150	7	60	5	0.035	5.5-6.5
EVD35-07-080-080	EVD35-07-080-120	EVD35-07-080-150	7	80	5	0.035	5.5-6.5
EVD35-07-100-080	EVD35-07-100-120	EVD35-07-100-150	7	100	5	0.035	5.5-6.5
EVD35-07-120-080	EVD35-07-120-120	EVD35-07-120-150	7	120	5	0.035	5.5-6.5
EVD35-07-150-080	EVD35-07-150-120	EVD35-07-150-150	7	150	5	0.035	5.5-6.5
EVD35-08-020-080	EVD35-08-020-120	EVD35-08-020-150	8	20	5	0.035	6.5-7.5
EVD35-08-040-080	EVD35-08-040-120	EVD35-08-040-150	8	40	5	0.035	6.5-7.5
EVD35-08-060-080	EVD35-08-060-120	EVD35-08-060-150	8	60	5	0.035	6.5-7.5
EVD35-08-080-080	EVD35-08-080-120	EVD35-08-080-150	8	80	5	0.035	6.5-7.5
EVD35-08-100-080	EVD35-08-100-120	EVD35-08-100-150	8	100	5	0.035	6.5-7.5
EVD35-08-120-080	EVD35-08-120-120	EVD35-08-120-150	8	120	5	0.035	6.5-7.5
EVD35-08-150-080	EVD35-08-150-120	EVD35-08-150-150	8	150	5	0.035	6.5-7.5

Each system includes one stent and delivery catheter system.

EverFlex™

Self-expanding Peripheral Stent System

When you need to stent, trust the precision, strength, and flexibility of the EverFlex stent,¹ indicated for use in the superficial femoral and proximal popliteal arteries (SFA/PPA), common iliac, and/or external iliac arteries. The EverFlex stent is available in lengths from 20 to 200 mm (SFA/PPA) or 20–120 mm (iliac). The broad size matrix, all deliverable through a 6 F catheter, provides the most appropriate single-stent fit. Peak-to-peak connection nodes disperse force uniformly, enhancing durability, while the spiral-cell connection and three-wave peak design optimize flexibility.¹

¹Data on file at Medtronic in reports R-2380, RE-PV10977, R-2398, and R-3348.

Catheter length		Stent size		Compatibility		
Reference number 80 cm	Reference number 120 cm	Diameter (mm)	Length (mm)	Recomm. introducer sheath (F)	Guidewire (in)	Recomm. lumen size (mm)
PRB35-06-020-080	PRB35-06-020-120	6	20	6	0.035	4.5–5.5
PRB35-06-030-080	PRB35-06-030-120	6	30	6	0.035	4.5–5.5
PRB35-06-040-080	PRB35-06-040-120	6	40	6	0.035	4.5–5.5
PRB35-06-060-080	PRB35-06-060-120	6	60	6	0.035	4.5–5.5
PRB35-06-080-080	PRB35-06-080-120	6	80	6	0.035	4.5–5.5
PRB35-06-100-080	PRB35-06-100-120	6	100	6	0.035	4.5–5.5
PRB35-06-120-080	PRB35-06-120-120	6	120	6	0.035	4.5–5.5
PRB35-06-150-080	PRB35-06-150-120	6	150	6	0.035	4.5–5.5
N/A	PRB35-06-200-120	6	200	6	0.035	4.5–5.5
PRB35-07-020-080	PRB35-07-020-120	7	20	6	0.035	5.5–6.5
PRB35-07-030-080	PRB35-07-030-120	7	30	6	0.035	5.5–6.5
PRB35-07-040-080	PRB35-07-040-120	7	40	6	0.035	5.5–6.5
PRB35-07-060-080	PRB35-07-060-120	7	60	6	0.035	5.5–6.5
PRB35-07-080-080	PRB35-07-080-120	7	80	6	0.035	5.5–6.5
PRB35-07-100-080	PRB35-07-100-120	7	100	6	0.035	5.5–6.5
PRB35-07-120-080	PRB35-07-120-120	7	120	6	0.035	5.5–6.5
PRB35-07-150-080	PRB35-07-150-120	7	150	6	0.035	5.5–6.5
N/A	PRB35-07-200-120	7	200	6	0.035	5.5–6.5
PRB35-08-020-080	PRB35-08-020-120	8	20	6	0.035	6.5–7.5
PRB35-08-030-080	PRB35-08-030-120	8	30	6	0.035	6.5–7.5
PRB35-08-040-080	PRB35-08-040-120	8	40	6	0.035	6.5–7.5
PRB35-08-060-080	PRB35-08-060-120	8	60	6	0.035	6.5–7.5
PRB35-08-080-080	PRB35-08-080-120	8	80	6	0.035	6.5–7.5
PRB35-08-100-080	PRB35-08-100-120	8	100	6	0.035	6.5–7.5
PRB35-08-120-080	PRB35-08-120-120	8	120	6	0.035	6.5–7.5
PRB35-08-150-080	PRB35-08-150-120	8	150	6	0.035	6.5–7.5
N/A	PRB35-08-200-120	8	200	6	0.035	6.5–7.5

Each system includes one stent and delivery catheter system.

Protégé™ EverFlex™

Self-expanding Biliary Stent System

The Protégé EverFlex stent system is designed for the palliative treatment of malignant neoplasms in the biliary tree. The stent is made of nitinol and comes pre-mounted on an over-the-wire delivery system. The proximal and distal ends of the stent have tantalum radiopaque markers for enhanced visibility.

Reference number	Unconstrained stent diameter (mm)	Unconstrained stent length (mm)	Rec. lumen size (mm)	Usable length (cm)	Sheath size (F)	Guidewire acceptance (in)	Outside diameter (in)
PRB35-05-020-080	5	20	3.5-4.5	80	6	0.035	0.079
PRB35-05-020-120	5	20	3.5-4.5	120	6	0.035	0.079
PRB35-05-030-080	5	30	3.5-4.5	80	6	0.035	0.079
PRB35-05-030-120	5	30	3.5-4.5	120	6	0.035	0.079
PRB35-05-040-080	5	40	3.5-4.5	80	6	0.035	0.079
PRB35-05-040-120	5	40	3.5-4.5	120	6	0.035	0.079
PRB35-05-060-080	5	60	3.5-4.5	80	6	0.035	0.079
PRB35-05-060-120	5	60	3.5-4.5	120	6	0.035	0.079
PRB35-05-080-080	5	80	3.5-4.5	80	6	0.035	0.079
PRB35-05-080-120	5	80	3.5-4.5	120	6	0.035	0.079
PRB35-05-100-080	5	100	3.5-4.5	80	6	0.035	0.079
PRB35-05-100-120	5	100	3.5-4.5	120	6	0.035	0.079
PRB35-05-120-080	5	120	3.5-4.5	80	6	0.035	0.079
PRB35-05-120-120	5	120	3.5-4.5	120	6	0.035	0.079

Each system includes one stent and delivery catheter system.



Protégé™ GPS™

Self-expanding Peripheral and Biliary Stent System

The Protégé GPS self-expanding peripheral and biliary stent system is a self-expanding, nitinol stent system indicated for the common and external iliac arteries (excluding 14 mm sizes), and designed for the palliative treatment of malignant neoplasms in the biliary tree. The stent is made of a nickel titanium alloy (nitinol) and comes pre-mounted on a 6 F over-the-wire delivery system. The stent is cut from a nitinol tube into an open lattice design and has tantalum radiopaque markers at the proximal and distal ends of the stent.

Reference number	Unconstrained stent diameter (mm)	Unconstrained stent length (mm)	Rec. lumen size (mm)	Usable length (cm)	Sheath size (F)	Guidewire acceptance (in)	Outside diameter (in)
SERB65-09-20-80	9	20	7.5-8.5	80	6	0.035	0.079
SERB65-09-30-80	9	30	7.5-8.5	80	6	0.035	0.079
SERB65-09-40-80	9	40	7.5-8.5	80	6	0.035	0.079
SERB65-09-60-80	9	60	7.5-8.5	80	6	0.035	0.079
SERB65-09-80-80	9	80	7.5-8.5	80	6	0.035	0.079
SERB65-10-20-80	10	20	8.5-9.5	80	6	0.035	0.079
SERB65-10-30-80	10	30	8.5-9.5	80	6	0.035	0.079
SERB65-10-40-80	10	40	8.5-9.5	80	6	0.035	0.079
SERB65-10-60-80	10	60	8.5-9.5	80	6	0.035	0.079
SERB65-10-80-80	10	80	8.5-9.5	80	6	0.035	0.079
SERB65-12-20-80	12	20	9.5-11.0	80	6	0.035	0.079
SERB65-12-30-80	12	30	9.5-11.0	80	6	0.035	0.079
SERB65-12-40-80	12	40	9.5-11.0	80	6	0.035	0.079
SERB65-12-60-80	12	60	9.5-11.0	80	6	0.035	0.079
SERB65-12-80-80	12	80	9.5-11.0	80	6	0.035	0.079
SERB65-09-20-120	9	20	7.5-8.5	120	6	0.035	0.079
SERB65-09-30-120	9	30	7.5-8.5	120	6	0.035	0.079
SERB65-09-40-120	9	40	7.5-8.5	120	6	0.035	0.079
SERB65-09-60-120	9	60	7.5-8.5	120	6	0.035	0.079
SERB65-09-80-120	9	80	7.5-8.5	120	6	0.035	0.079
SERB65-10-20-120	10	20	8.5-9.5	120	6	0.035	0.079
SERB65-10-30-120	10	30	8.5-9.5	120	6	0.035	0.079
SERB65-10-40-120	10	40	8.5-9.5	120	6	0.035	0.079
SERB65-10-60-120	10	60	8.5-9.5	120	6	0.035	0.079
SERB65-10-80-120	10	80	8.5-9.5	120	6	0.035	0.079
SERB65-12-20-120	12	20	9.5-11.0	120	6	0.035	0.079
SERB65-12-30-120	12	30	9.5-11.0	120	6	0.035	0.079
SERB65-12-40-120	12	40	9.5-11.0	120	6	0.035	0.079
SERB65-12-60-120	12	60	9.5-11.0	120	6	0.035	0.079
SERB65-12-80-120	12	80	9.5-11.0	120	6	0.035	0.079

Biliary Only

SERB65-14-20-80	14	20	11.5-13.0	80	6	0.035	0.079
SERB65-14-30-80	14	30	11.5-13.0	80	6	0.035	0.079
SERB65-14-40-80	14	40	11.5-13.0	80	6	0.035	0.079
SERB65-14-60-80	14	60	11.5-13.0	80	6	0.035	0.079
SERB65-14-80-80	14	80	11.5-13.0	80	6	0.035	0.079
SERB65-14-20-120	14	20	11.5-13.0	120	6	0.035	0.079
SERB65-14-30-120	14	30	11.5-13.0	120	6	0.035	0.079
SERB65-14-40-120	14	40	11.5-13.0	120	6	0.035	0.079
SERB65-14-60-120	14	60	11.5-13.0	120	6	0.035	0.079
SERB65-14-80-120	14	80	11.5-13.0	120	6	0.035	0.079

Each system includes one stent and delivery catheter system.

Protégé™ RX

Self-expanding Carotid Stent System

The Protégé RX stent system is designed for carotid artery stenting. The nitinol stent comes pre-mounted on a 6 F, 0.014" rapid exchange delivery system. The proximal and distal ends of the stent have tantalum radiopaque markers for enhanced visibility.

Reference number	Unconstrained stent diameter (mm)	Unconstrained stent length (mm)	Recommended lumen size (mm)	Usable catheter length (cm)	Sheath size (F)	Guidewire acceptance (in)
Straight						
SECX-6-20-135	6	20	4.5-5.5	135	6	0.014
SECX-6-30-135	6	30	4.5-5.5	135	6	0.014
SECX-6-40-135	6	40	4.5-5.5	135	6	0.014
SECX-6-60-135	6	60	4.5-5.5	135	6	0.014
SECX-7-20-135	7	20	5.5-6.5	135	6	0.014
SECX-7-30-135	7	30	5.5-6.5	135	6	0.014
SECX-7-40-135	7	40	5.5-6.5	135	6	0.014
SECX-7-60-135	7	60	5.5-6.5	135	6	0.014
SECX-8-20-135	8	20	6.5-7.5	135	6	0.014
SECX-8-30-135	8	30	6.5-7.5	135	6	0.014
SECX-8-40-135	8	40	6.5-7.5	135	6	0.014
SECX-8-60-135	8	60	6.5-7.5	135	6	0.014
SECX-9-20-135	9	20	7.5-8.5	135	6	0.014
SECX-9-30-135	9	30	7.5-8.5	135	6	0.014
SECX-9-40-135	9	40	7.5-8.5	135	6	0.014
SECX-9-60-135	9	60	7.5-8.5	135	6	0.014
SECX-10-20-135	10	20	8.5-9.5	135	6	0.014
SECX-10-30-135	10	30	8.5-9.5	135	6	0.014
SECX-10-40-135	10	40	8.5-9.5	135	6	0.014
SECX-10-60-135	10	60	8.5-9.5	135	6	0.014
Tapered						
SECX-8-6-30-135	8/6	30	(6.5-7.5)-(4.5-5.5)	135	6	0.014
SECX-8-6-40-135	8/6	40	(6.5-7.5)-(4.5-5.5)	135	6	0.014
SECX-10-7-30-135	10/7	30	(8.5-9.5)-(5.5-6.5)	135	6	0.014
SECX-10-7-40-135	10/7	40	(8.5-9.5)-(5.5-6.5)	135	6	0.014

Each system includes one stent and delivery catheter system.

Visi-Pro™

Balloon-expandable Peripheral Stent System

The Visi-Pro stent system is indicated for use in the common and external iliac, and biliary arteries. It is made from a stainless steel tube that is cut into an open lattice design and has tantalum radiopaque markers at the proximal and distal ends of the stent. The stent is premounted on a balloon catheter delivery system, and is expanded and deployed by inflating the balloon.

Catheter length		Stent size		Balloon	Compatibility	
Reference number 80 cm	Reference number 135 cm	Diameter (mm)	Length (mm)	Balloon length (mm)	Rec. introducer sheath (F)	Guidewire (in)
PXB35-05-12-080	-	5	12	15	6†	0.035
PXB35-05-17-080	PXB35-05-17-135	5	17	20	6†	0.035
PXB35-05-27-080	PXB35-05-27-135	5	27	30	6†	0.035
PXB35-05-37-080	PXB35-05-37-135	5	37	40	6†	0.035
PXB35-05-57-080	PXB35-05-57-135	5	57	60	6†	0.035
PXB35-06-12-080	-	6	12	15	6†	0.035
PXB35-06-17-080	PXB35-06-17-135	6	17	20	6†	0.035
PXB35-06-27-080	PXB35-06-27-135	6	27	30	6†	0.035
PXB35-06-37-080	PXB35-06-37-135	6	37	40	6†	0.035
PXB35-06-57-080	PXB35-06-57-135	6	57	60	6†	0.035
PXB35-07-12-080	-	7	12	15	6†	0.035
PXB35-07-17-080	PXB35-07-17-135	7	17	20	6†	0.035
PXB35-07-27-080	PXB35-07-27-135	7	27	30	6†	0.035
PXB35-07-37-080	PXB35-07-37-135	7	37	40	6†	0.035
PXB35-07-57-080	PXB35-07-57-135	7	57	60	6†	0.035
PXB35-08-17-080	PXB35-08-17-135	8	17	20	6†	0.035
PXB35-08-27-080	PXB35-08-27-135	8	27	30	6†	0.035
PXB35-08-37-080	PXB35-08-37-135	8	37	40	6†	0.035
PXB35-08-57-080	PXB35-08-57-135	8	57	60	6†	0.035
PXB35-09-17-080	PXB35-09-17-135	9	17	20	7	0.035
PXB35-09-27-080	PXB35-09-27-135	9	27	30	7	0.035
PXB35-09-37-080	PXB35-09-37-135	9	37	40	7	0.035
PXB35-09-57-080	PXB35-09-57-135	9	57	60	7	0.035
PXB35-10-17-080	PXB35-10-17-135	10	17	20	7	0.035
PXB35-10-27-080	PXB35-10-27-135	10	27	30	7	0.035
PXB35-10-37-080	PXB35-10-37-135	10	37	40	7	0.035
PXB35-10-57-080	PXB35-10-57-135	10	57	60	7	0.035

†6 F = 0.085" I.D

Each system includes one stent and delivery catheter system.

IntraStent™

Biliary Stent System

The IntraStent biliary stent is made from a stainless steel tube that is cut into an open lattice design. Its design allows it to be crimped onto a non-compliant PTA balloon catheter. After it is mounted onto a balloon catheter, the stent is expanded and deployed by inflating the balloon.

Reference number	Unexpanded stent size		Expanded stent size	
	Diameter (mm)	Length (mm)	Diameter (mm)	Length (mm)
IntraStent™ DoubleStrut LD Biliary Stent				
S15-16 (90-1431-000)	3.8	16.0	9, 10, 11, 12	16.0
S15-26 (90-1431-001)	3.8	26.0	9, 10, 11, 12	26.0
S15-36 (90-1431-002)	3.8	36.0	9, 10, 11, 12	36.0
S15-56 (90-1431-003)	3.8	56.0	9, 10, 11, 12	56.0
S15-76 (90-1431-004)	3.8	76.0	9, 10, 11, 12	76.0
IntraStent™ Mega LD Biliary Stent				
S17-16 (90-2313-000)	3.8	16.0	9, 10, 12	16.0
S17-26 (90-2313-001)	3.8	26.0	9, 10, 12	26.0
S17-36 (90-2313-002)	3.8	36.0	9, 10, 12	36.0
IntraStent™ Max LD Biliary Stent				
S18-16 (90-2319-000)	4.5	16.0	12	16.0
S18-26 (90-2319-001)	4.5	26.0	12	26.0
S18-36 (90-2319-002)	4.5	36.0	12	36.0

Abre™

Venous Self-expanding Stent System

The Abre venous self-expanding stent system is designed for the unique challenges of venous disease. It offers easy deployment to let physicians focus on their patient, and delivers demonstrated endurance to give patients freedom of movement.^{1,2}

¹ Abre stent IFU.

² Test data on file at Medtronic. Report 10558227DOC_Rev A. Bench test results may not be indicative of clinical performance.

Stent diameters	Stent Lengths and Product Numbers					
	40 mm	60 mm	80 mm	100 mm	120 mm	150 mm
10 mm	AB9U10040090	AB9U10060090	AB9U10080090	AB9U10100090	AB9U10120090	AB9U10150090
12 mm	-	AB9U12060090	AB9U12080090	AB9U12100090	AB9U12120090	AB9U12150090
14 mm	-	AB9U14060090	AB9U14080090	AB9U14100090	AB9U14120090	AB9U14150090
16 mm	-	AB9U16060090	AB9U16080090	AB9U16100090	AB9U16120090	AB9U16150090
18 mm	-	AB9U18060090	AB9U18080090	AB9U18100090	AB9U18120090	AB9U18150090
20 mm	-	AB9U20060090	AB9U20080090	AB9U20100090	AB9U20120090	AB9U20150090

Risks of the Abre system may include pain, myocardial infarction, pulmonary embolism, and restenosis of stented segment.

PTA Balloons



Fortrex™
HP PTA Balloon Catheter



Pacific™ Plus
PTA Catheter

0.035" OTW PTA Dilatation Catheter

The 0.035" EverCross balloon catheter is available in a broad range of balloon sizes.

Reference number usable length 40 cm	Reference number usable length 80 cm	Reference number usable length 135 cm	Balloon diameter (mm)	Balloon length (mm)	Recomm. introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
-	AB35W03020080	AB35W03020135	3.0	20	5	10	20
-	AB35W03030080	AB35W03030135	3.0	30	5	10	20
-	AB35W03040080	AB35W03040135	3.0	40	5	10	20
-	AB35W03060080	AB35W03060135	3.0	60	5	10	20
-	AB35W03080080	AB35W03080135	3.0	80	5	10	20
-	AB35W03100080	AB35W03100135	3.0	100	5	10	20
-	AB35W03120080	AB35W03120135	3.0	120	5	10	20
-	AB35W03150080	AB35W03150135	3.0	150	5	10	20
-	AB35W03200080	AB35W03200135	3.0	200	5	10	20
-	AB35W04020080	AB35W04020135	4.0	20	5	10	20
-	AB35W04030080	AB35W04030135	4.0	30	5	10	20
-	AB35W04040080	AB35W04040135	4.0	40	5	10	20
-	AB35W04060080	AB35W04060135	4.0	60	5	10	20
-	AB35W04080080	AB35W04080135	4.0	80	5	10	20
-	AB35W04100080	AB35W04100135	4.0	100	5	10	20
-	AB35W04120080	AB35W04120135	4.0	120	5	10	20
-	AB35W04150080	AB35W04150135	4.0	150	5	10	20
-	AB35W04200080	AB35W04200135	4.0	200	5	10	20
AB35W05020040	AB35W05020080	AB35W05020135	5.0	20	5	10	18
AB35W05030040	AB35W05030080	AB35W05030135	5.0	30	5	10	18
AB35W05040040	AB35W05040080	AB35W05040135	5.0	40	5	10	18
AB35W05060040	AB35W05060080	AB35W05060135	5.0	60	5	10	18
AB35W05080040	AB35W05080080	AB35W05080135	5.0	80	5	10	18
-	AB35W05100080	AB35W05100135	5.0	100	5	10	18
AB35W05120040	AB35W05120080	AB35W05120135	5.0	120	5	10	16
-	AB35W05150080	AB35W05150135	5.0	150	5	10	16
-	AB35W05200080	AB35W05200135	5.0	200	5	10	16
AB35W06020040	AB35W06020080	AB35W06020135	6.0	20	5	8	14
-	AB35W06030080	AB35W06030135	6.0	30	5	8	14
AB35W06040040	AB35W06040080	AB35W06040135	6.0	40	5	8	14
-	AB35W06060080	AB35W06060135	6.0	60	5	8	14
AB35W06080040	AB35W06080080	AB35W06080135	6.0	80	5	8	14
-	AB35W06100080	AB35W06100135	6.0	100	5	8	14
AB35W06120040	AB35W06120080	AB35W06120135	6.0	120	5	8	12
-	AB35W06150080	AB35W06150135	6.0	150	5	8	12
-	AB35W06200080	AB35W06200135	6.0	200	6	8	11
AB35W07020040	AB35W07020080	AB35W07020135	7.0	20	5	7	14
-	AB35W07030080	AB35W07030135	7.0	30	5	7	14
AB35W07040040	AB35W07040080	AB35W07040135	7.0	40	5	7	14
AB35W07060040	AB35W07060080	AB35W07060135	7.0	60	6	7	14
-	AB35W07080080	AB35W07080135	7.0	80	6	7	14
-	AB35W07100080	AB35W07100135	7.0	100	6	7	14
-	AB35W07120080	AB35W07120135	7.0	120	6	7	10
-	AB35W07150080	AB35W07150135	7.0	150	6	7	10
-	AB35W07200080	AB35W07200135	7.0	200	6	7	10
AB35W08020040	AB35W08020080	AB35W08020135	8.0	20	6	7	14
-	AB35W08030080	AB35W08030135	8.0	30	6	7	14
AB35W08040040	AB35W08040080	AB35W08040135	8.0	40	6	7	14
AB35W08060040	AB35W08060080	AB35W08060135	8.0	60	6	7	14
-	AB35W08080080	AB35W08080135	8.0	80	6	7	14
-	AB35W09020080	AB35W09020135	9.0	20	6	7	12
-	AB35W09030080	AB35W09030135	9.0	30	6	7	12
-	AB35W09040080	AB35W09040135	9.0	40	6	7	12
-	AB35W09060080	AB35W09060135	9.0	60	6	7	12
-	AB35W09080080	AB35W09080135	9.0	80	6	7	12
-	AB35W10020080	AB35W10020135	10.0	20	6	7	11
-	AB35W10030080	AB35W10030135	10.0	30	6	7	11
-	AB35W10040080	AB35W10040135	10.0	40	6	7	11
-	AB35W10060080	AB35W10060135	10.0	60	7	7	11
-	AB35W12020080	AB35W12020135	12.0	20	7	7	10
-	AB35W12040080	AB35W12040135	12.0	40	7	7	10
-	AB35W12060080	AB35W12060135	12.0	60	7	7	10

Admiral Xtreme™

PTA Balloon Catheter

Treat longer peripheral arterial disease (PAD) lesions above the knee with the Admiral Xtreme PTA balloon catheter. The Admiral Xtreme is offered in 250 and 300 mm lengths, with 0.035" guidewire compatibility, featuring large inflation lumen for rapid inflation and deflation.†

†Data on file at Medtronic in 10613058DOC, CP076TP02R, and CP076TP04R. Bench test results may not be indicative of clinical performance.

Reference number usable length 80 cm	Reference number usable length 130 cm	Balloon diameter (mm)	Balloon length (mm)	Recom. introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
ADM040250080	ADM040250130	4.0	250	5	6	14
ADM040300080	ADM040300130	4.0	300	5	6	14
ADM050250080	ADM050250130	5.0	250	5	6	14
ADM050300080L	ADM050300130L	5.0	300	5	6	14
ADM060250080L	ADM060250130L	6.0	250	5	6	12
ADM060300080L	ADM060300130L	6.0	300	5	6	12
ADM070250080	ADM070250130	7.0	250	6	6	12

Fortrex™

0.035" OTW HP PTA Balloon Catheter

The Fortrex PTA balloon is the high pressure solution for AV access, and it is also intended for use in the peripheral vascular system. Engineered specifically for peak performance at rated burst pressure, the Fortrex balloon offers the deliverability, predictability, and procedural efficiency† desired by physicians.

†Test data on file at Medtronic in report RE-PV14089. Bench test results may not be indicative of clinical performance.

Reference number usable length 40 cm	Reference number usable length 80 cm	Reference number usable length 135 mm	Balloon diameter (mm)	Balloon length (mm)	Nominal pressure (atm)	Rated burst pressure (atm)	Recommended introducer sheath (F)
A35HPV04020040	A35HPV04020080	A35HPV04020135	4.0	20	12	24	6
A35HPV04040040	A35HPV04040080	A35HPV04040135	4.0	40	12	24	6
A35HPV04080040	A35HPV04080080	A35HPV04080135	4.0	80	12	24	6
A35HPV04100040	A35HPV04100080	A35HPV04100135	4.0	100	12	24	6
A35HPV05020040	A35HPV05020080	A35HPV05020135	5.0	20	12	24	6
A35HPV05040040	A35HPV05040080	A35HPV05040135	5.0	40	12	24	6
A35HPV05080040	A35HPV05080080	A35HPV05080135	5.0	80	12	24	6
A35HPV05100040	A35HPV05100080	A35HPV05100135	5.0	100	12	24	6
A35HPV06020040	A35HPV06020080	A35HPV06020135	6.0	20	12	24	6
A35HPV06040040	A35HPV06040080	A35HPV06040135	6.0	40	12	24	6
A35HPV06080040	A35HPV06080080	A35HPV06080135	6.0	80	12	23	6
A35HPV06100040	A35HPV06100080	A35HPV06100135	6.0	100	12	23	6
A35HPV07020040	A35HPV07020080	A35HPV07020135	7.0	20	9	20	6
A35HPV07040040	A35HPV07040080	A35HPV07040135	7.0	40	9	20	6
A35HPV07080040	A35HPV07080080	A35HPV07080135	7.0	80	9	20	6
A35HPV07100040	A35HPV07100080	A35HPV07100135	7.0	100	9	20	6
A35HPV08040040	A35HPV08040080	A35HPV08040135	8.0	40	9	20	6
A35HPV08080040	A35HPV08080080	A35HPV08080135	8.0	80	9	19	6
A35HPV08100040	A35HPV08100080	A35HPV08100135	8.0	100	9	18	6
A35HPV09040040	A35HPV09040080	A35HPV09040135	9.0	40	9	18	7
A35HPV09080040	A35HPV09080080	A35HPV09080135	9.0	80	9	18	7
A35HPV10040040	A35HPV10040080	A35HPV10040135	10.0	40	8	16	7
A35HPV10080040	A35HPV10080080	A35HPV10080135	10.0	80	8	16	7
A35HPV12040040	A35HPV12040080	A35HPV12040135	12.0	40	8	14	7
A35HPV12080040	A35HPV12080080	A35HPV12080135	12.0	80	7	12	7

Risks of Fortrex balloon may include abrupt or subacute closure, pseudo-aneurysm, thrombosis, allergic reaction to device materials or procedural medications, arterial injury (such as dissection, perforation or rupture), embolism, hematoma, and hemorrhage.

Pacific™ Plus

PTA Catheter OTW 0.018"



The Pacific Plus PTA catheter is the go-to balloon for treating peripheral arterial disease (PAD) lesions. Through its broad size matrix and low 0.018" profile, it offers versatility to handle cases above the knee (ATK) and below the knee (BTK) from various access points.

Usable length 90 cm	Usable length 130 cm	Usable length 150 cm	Usable length 180 cm	Usable length 200 cm	Balloon diameter (mm)	Balloon length (mm)	Minimum sheath inner diameter (F)	Nominal pressure (atm)	RBP (atm)
PCP020020090	PCP020020130	PCP020020150	PCP020020180	-	2.00	20	4	8	22
PCP020040090	PCP020040130	PCP020040150	PCP020040180	-	2.00	40	4	8	22
PCP020060090	PCP020060130	PCP020060150	PCP020060180	-	2.00	60	4	8	22
PCP020080090	PCP020080130	PCP020080150	PCP020080180	-	2.00	80	4	8	22
PCP020120090	PCP020120130	PCP020120150	PCP020120180	-	2.00	120	4	8	22
PCP020150090	PCP020150130	PCP020150150	PCP020150180	-	2.00	150	4	8	22
PCP025020090	PCP025020130	PCP025020150	PCP025020180	-	2.50	20	4	8	16
PCP025040090	PCP025040130	PCP025040150	PCP025040180	-	2.50	40	4	8	16
PCP025060090	PCP025060130	PCP025060150	PCP025060180	-	2.50	60	4	8	16
PCP025080090	PCP025080130	PCP025080150	PCP025080180	-	2.50	80	4	8	16
PCP025120090	PCP025120130	PCP025120150	PCP025120180	-	2.50	120	4	8	16
PCP025150090	PCP025150130	PCP025150150	PCP025150180	-	2.50	150	4	8	16
PCP030020090	PCP030020130	PCP030020150	PCP030020180	-	3.00	20	4	8	16
PCP030040090	PCP030040130	PCP030040150	PCP030040180	-	3.00	40	4	8	16
PCP030060090	PCP030060130	PCP030060150	PCP030060180	-	3.00	60	4	8	16
PCP030080090	PCP030080130	PCP030080150	PCP030080180	-	3.00	80	4	8	16
PCP030120090	PCP030120130	PCP030120150	PCP030120180	-	3.00	120	4	8	16
PCP030150090	PCP030150130	PCP030150150	PCP030150180	-	3.00	150	4	8	16
PCP035020090	PCP035020130	PCP035020150	PCP035020180	-	3.50	20	4	8	16
PCP035040090	PCP035040130	PCP035040150	PCP035040180	-	3.50	40	4	8	16
PCP035060090	PCP035060130	PCP035060150	PCP035060180	-	3.50	60	4	8	16
PCP035080090	PCP035080130	PCP035080150	PCP035080180	-	3.50	80	4	8	16
PCP035120090	PCP035120130	PCP035120150	PCP035120180	-	3.50	120	4	8	16
PCP035150090	PCP035150130	PCP035150150	PCP035150180	-	3.50	150	4	8	16
PCP040020090	-	PCP040020150	-	PCP040020200	4.00	20	4	8	14
PCP040040090	-	PCP040040150	-	PCP040040200	4.00	40	4	8	14
PCP040060090	-	PCP040060150	-	PCP040060200	4.00	60	4	8	14
PCP040080090	-	PCP040080150	-	PCP040080200	4.00	80	4	8	14
PCP040100090	-	PCP040100150	-	PCP040100200	4.00	100	4	8	14
PCP040120090	-	PCP040120150	-	PCP040120200	4.00	120	4	8	14
PCP040150090	-	PCP040150150	-	PCP040150200	4.00	150	4	8	14
PCP040200090	-	PCP040200150	-	PCP040200200	4.00	200	4	8	14
PCP040250090	-	PCP040250150	-	PCP040250200	4.00	250	4	8	14
PCP040300090	-	PCP040300150	-	PCP040300200	4.00	300	4	8	14
PCP050020090	-	PCP050020150	-	PCP050020200	5.00	20	4	8	14
PCP050040090	-	PCP050040150	-	PCP050040200	5.00	40	4	8	14
PCP050060090	-	PCP050060150	-	PCP050060200	5.00	60	4	8	14
PCP050080090	-	PCP050080150	-	PCP050080200	5.00	80	4	8	14
PCP050100090	-	PCP050100150	-	PCP050100200	5.00	100	4	8	14
PCP050120090	-	PCP050120150	-	PCP050120200	5.00	120	4	8	14
PCP050150090	-	PCP050150150	-	PCP050150200	5.00	150	4	8	14
PCP050200090	-	PCP050200150	-	PCP050200200	5.00	200	4	8	14

Usable length 90 cm	Usable length 130 cm	Usable length 150 cm	Usable length 180 cm	Usable length 200 cm	Balloon diameter (mm)	Balloon length (mm)	Minimum sheath inner diameter (F)	Nominal pressure (atm)	RBP (atm)
PCP050300090	-	PCP050300150	-	PCP050300200	5.00	300	5	8	14
PCP060020090	-	PCP060020150	-	PCP060020200	6.00	20	4	8	14
PCP060040090	-	PCP060040150	-	PCP060040200	6.00	40	4	8	14
PCP060060090	-	PCP060060150	-	PCP060060200	6.00	60	4	8	14
PCP060080090	-	PCP060080150	-	PCP060080200	6.00	80	4	8	14
PCP060100090	-	PCP060100150	-	PCP060100200	6.00	100	4	8	14
PCP060120090	-	PCP060120150	-	PCP060120200	6.00	120	4	8	14
PCP060150090	-	PCP060150150	-	PCP060150200	6.00	150	5	8	12
PCP060200090	-	PCP060200150	-	PCP060200200	6.00	200	5	8	12
PCP060250090	-	PCP060250150	-	PCP060250200	6.00	250	5	8	12
PCP060300090	-	PCP060300150	-	PCP060300200	6.00	300	5	8	12
PCP070020090	-	PCP070020150	-	PCP070020200	7.00	20	4	8	12
PCP070040090	-	PCP070040150	-	PCP070040200	7.00	40	4	8	12
PCP070060090	-	PCP070060150	-	PCP070060200	7.00	60	4	8	12
PCP070080090	-	PCP070080150	-	PCP070080200	7.00	80	5	8	12
PCP070100090	-	PCP070100150	-	PCP070100200	7.00	100	5	8	12
PCP070120090	-	PCP070120150	-	PCP070120200	7.00	120	5	8	12
PCP070150090	-	PCP070150150	-	PCP070150200	7.00	150	5	8	12
PCP070200090	-	PCP070200150	-	PCP070200200	7.00	200	5	8	12
PCP070250090	-	PCP070250150	-	PCP070250200	7.00	250	5	8	12

NanoCross™ Elite

0.014" OTW PTA Balloon Dilatation Catheter

Go deep and cross difficult lesions with NanoCross Elite balloon catheter.[†] Its seamless design, from balloon to tip, offers efficient energy transfer.[‡]

[†]Medtronic data on file in reports RE-PV12779, RE-PV13497.

[‡]Medtronic data on file in reports RE-PV12779, RE-PV13497, RE-PV13040, RE-PV13834, RE-PV12882, and RE-PV13697.

Reference number usable length 90 cm	Reference number usable length 150 cm	Balloon diameter (mm)	Balloon length (mm)	Recomm. introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
AB14W015020090	AB14W015020150	1.5	20	4	8	14
AB14W015040090	AB14W015040150	1.5	40	4	8	14
AB14W020020090	AB14W020020150	2.0	20	4	8	14
AB14W020040090	AB14W020040150	2.0	40	4	8	14
AB14W020060090	AB14W020060150	2.0	60	4	8	14
AB14W020080090	AB14W020080150	2.0	80	4	8	14
AB14W020100090	AB14W020100150	2.0	100	4	8	14
AB14W020120090	AB14W020120150	2.0	120	4	8	14
AB14W020150090	AB14W020150150	2.0	150	4	8	14
AB14W020210090	AB14W020210150	2.0 (proximal)/1.5 (distal)	210	4	8	14
AB14W025020090	AB14W025020150	2.5	20	4	8	14
AB14W025040090	AB14W025040150	2.5	40	4	8	14
AB14W025060090	AB14W025060150	2.5	60	4	8	14
AB14W025080090	AB14W025080150	2.5	80	4	8	14
AB14W025100090	AB14W025100150	2.5	100	4	8	14
AB14W025120090	AB14W025120150	2.5	120	4	8	14
AB14W025150090	AB14W025150150	2.5	150	4	8	14
AB14W025210090	AB14W025210150	2.5 (proximal)/2.0 (distal)	210	4	8	14
AB14W030020090	AB14W030020150	3.0	20	4	8	14
AB14W030040090	AB14W030040150	3.0	40	4	8	14
AB14W030060090	AB14W030060150	3.0	60	4	8	14
AB14W030080090	AB14W030080150	3.0	80	4	8	14
AB14W030100090	AB14W030100150	3.0	100	4	8	14
AB14W030120090	AB14W030120150	3.0	120	4	8	14
AB14W030150090	AB14W030150150	3.0	150	4	8	14
AB14W030210090	AB14W030210150	3.0 (proximal)/2.5 (distal)	210	4	8	14
AB14W035020090	AB14W035020150	3.5	20	4	8	14
AB14W035040090	AB14W035040150	3.5	40	4	8	14
AB14W035060090	AB14W035060150	3.5	60	4	8	14
AB14W035080090	AB14W035080150	3.5	80	4	8	14
AB14W035100090	AB14W035100150	3.5	100	4	8	14
AB14W035120090	AB14W035120150	3.5	120	4	8	14
AB14W035150090	AB14W035150150	3.5	150	4	8	14
AB14W035210090	AB14W035210150	3.5 (proximal)/3.0 (distal)	210	4	8	14
AB14W040020090	AB14W040020150	4.0	20	4	8	14
AB14W040040090	AB14W040040150	4.0	40	4	8	14
AB14W040060090	AB14W040060150	4.0	60	4	8	14
AB14W040080090	AB14W040080150	4.0	80	4	8	14
AB14W040100090	AB14W040100150	4.0	100	4	8	14
AB14W040120090	AB14W040120150	4.0	120	4	8	14
AB14W040150090	AB14W040150150	4.0	150	4	8	14
AB14W040210090	AB14W040210150	4.0 (proximal)/3.5 (distal)	210	4	8	14
AB14W050020090	AB14W050020150	5.0	20	5	8	14
AB14W050040090	AB14W050040150	5.0	40	5	8	14
AB14W050060090	AB14W050060150	5.0	60	5	8	14

Risks of NanoCross Elite balloon may include abrupt or subacute closure, pseudo-aneurysm, thrombosis, allergic reaction to device materials or procedural medications, arterial injury (such as dissection, perforation or rupture), embolism, hematoma, and hemorrhage.

Reference number usable length 90 cm	Reference number usable length 150 cm	Balloon diameter (mm)	Balloon length (mm)	Recomm. introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
AB14W050080090	AB14W050080150	5.0	80	5	8	14
AB14W050100090	AB14W050100150	5.0	100	5	8	14
AB14W050120090	AB14W050120150	5.0	120	5	8	14
AB14W050150090	AB14W050150150	5.0	150	5	8	14
AB14W050200090	AB14W050200150	5.0	200	5	8	14
AB14W060020090	AB14W060020150	6.0	20	5	8	14
AB14W060040090	AB14W060040150	6.0	40	5	8	14
AB14W060060090	AB14W060060150	6.0	60	5	8	14
AB14W060080090	AB14W060080150	6.0	80	5	8	14
AB14W060100090	AB14W060100150	6.0	100	5	8	14
AB14W060120090	AB14W060120150	6.0	120	5	8	14
AB14W060150090	AB14W060150150	6.0	150	5	8	14
AB14W060200090	AB14W060200150	6.0	200	6	8	14

RapidCross™

PTA Rapid Exchange Balloon Dilatation Catheter

Choose the RapidCross PTA rapid exchange balloon dilatation catheter when you need to reach more distal lesions in the small vessels below the knee. This catheter is the Medtronic solution for a rapid exchange balloon with 0.014" guidewire capability.

Reference number usable length 90 cm	Reference number usable length 170 cm	Balloon diameter (mm)	Balloon length (mm)	Recomm. introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
A14BX020020090	A14BX020020170	2.0	20	4	8	14
A14BX020040090	A14BX020040170	2.0	40	4	8	14
A14BX020060090	A14BX020060170	2.0	60	4	8	14
A14BX020080090	A14BX020080170	2.0	80	4	8	14
A14BX020100090	A14BX020100170	2.0	100	4	8	14
A14BX020120090	A14BX020120170	2.0	120	4	8	14
A14BX020150090	A14BX020150170	2.0	150	4	8	14
A14BX020210090	A14BX020210170	2.0 (proximal)/1.5 (distal)	210	4	8	14
A14BX025020090	A14BX025020170	2.5	20	4	8	14
A14BX025040090	A14BX025040170	2.5	40	4	8	14
A14BX025060090	A14BX025060170	2.5	60	4	8	14
A14BX025080090	A14BX025080170	2.5	80	4	8	14
A14BX025100090	A14BX025100170	2.5	100	4	8	14
A14BX025120090	A14BX025120170	2.5	120	4	8	14
A14BX025150090	A14BX025150170	2.5	150	4	8	14
A14BX025210090	A14BX025210170	2.5 (proximal)/2.0 (distal)	210	4	8	14
A14BX030020090	A14BX030020170	3.0	20	4	8	14
A14BX030040090	A14BX030040170	3.0	40	4	8	14
A14BX030060090	A14BX030060170	3.0	60	4	8	14
A14BX030080090	A14BX030080170	3.0	80	4	8	14
A14BX030100090	A14BX030100170	3.0	100	4	8	14
A14BX030120090	A14BX030120170	3.0	120	4	8	14
A14BX030150090	A14BX030150170	3.0	150	4	8	14
A14BX030210090	A14BX030210170	3.0 (proximal)/2.5 (distal)	210	4	8	14
A14BX035020090	A14BX035020170	3.5	20	4	8	14
A14BX035040090	A14BX035040170	3.5	40	4	8	14
A14BX035060090	A14BX035060170	3.5	60	4	8	14
A14BX035080090	A14BX035080170	3.5	80	4	8	14
A14BX035100090	A14BX035100170	3.5	100	4	8	14
A14BX035120090	A14BX035120170	3.5	120	4	8	14
A14BX035150090	A14BX035150170	3.5	150	4	8	14
A14BX035210090	A14BX035210170	3.5 (proximal)/3.0 (distal)	210	4	8	14
A14BX040020090	A14BX040020170	4.0	20	4	8	14
A14BX040040090	A14BX040040170	4.0	40	4	8	14
A14BX040060090	A14BX040060170	4.0	60	4	8	14
A14BX040080090	A14BX040080170	4.0	80	4	8	14
A14BX040100090	A14BX040100170	4.0	100	4	8	14
A14BX040120090	A14BX040120170	4.0	120	4	8	14
A14BX040150090	A14BX040150170	4.0	150	4	8	14
A14BX040210090	A14BX040210170	4.0 (proximal)/3.5 (distal)	210	4	8	14

PTA Specialty Balloon



Chocolate™
PTA Balloon Catheter

Chocolate™*

PTA Balloon Catheter

The Chocolate™ PTA balloon is designed to minimize vessel trauma, dissections, and the need for bailout stenting above or below the knee.† The balloon's unique nitinol constraining structure creates pillows and grooves that provide a predictable, uniform, and atraumatic dilatation.


†Medtronic data on file: CLR782 Chocolate Bar Report.

Reference number	Balloon diameter (mm)	Balloon length (mm)	Catheter length (cm)	Guidewire (in)	Introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
CB1415025040OTW	2.5	40	150	0.014"	5	9	14
CB1415025080OTW	2.5	80	150	0.014"	5	9	14
CB1415025120OTW	2.5	120	150	0.014"	5	9	14
CB1415030040OTW	3.0	40	150	0.014"	5	9	14
CB1415030080OTW	3.0	80	150	0.014"	5	9	14
CB1415030120OTW	3.0	120	150	0.014"	5	9	14
CB1413535040OTW	3.5	40	135	0.014"	5	9	14
CB1413535080OTW	3.5	80	135	0.014"	5	9	14
CB1413535120OTW	3.5	120	135	0.014"	5	9	14
CB1413540040OTW	4.0	40	135	0.014"	5	9	14
CB1413540080OTW	4.0	80	135	0.014"	5	9	14
CB1413540120OTW	4.0	120	135	0.014"	5	9	14
CB1812050040OTW	5.0	40	120	0.018"	6	6	12
CB1812050080OTW	5.0	80	120	0.018"	6	6	12
CB1812050120OTW	5.0	120	120	0.018"	6	6	12
CB1812060040OTW	6.0	40	120	0.018"	6	6	12
CB1812060080OTW	6.0	80	120	0.018"	6	6	12
CB1812060120OTW	6.0	120	120	0.018"	6	6	12

™Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Risks of Chocolate balloon may include abrupt or subacute closure, pseudo-aneurysm, thrombosis, allergic reaction to device materials or procedural medications, arterial injury (such as dissection, perforation or rupture), embolism, hematoma, and hemorrhage.

Infusion Therapy



Cragg-McNamara™
Valved Infusion Catheters

Cragg-McNamara™

Valved Infusion Catheter

The Cragg-McNamara catheter is a 4 or 5 F single-lumen catheter with a proprietary valved tip, which gives the option to infuse without a guidewire in place.

Reference number	Introducer sheath (F)	Usable length (cm)	Infusion length (cm)	Max. guidewire (in)
41032-01	4	40	10	0.035
41033-01	4	40	20	0.035
41034-01	4	65	5	0.035
41035-01	4	65	10	0.035
41036-01	4	65	20	0.035
41037-01	4	100	5	0.035
41038-01	4	100	10	0.035
41039-01	4	100	20	0.035
41040-01	4	135	5	0.035
41041-01	4	135	10	0.035
41042-01	4	135	20	0.035
41043-01	5	40	5	0.038
41044-01	5	40	10	0.038
41045-01	5	40	20	0.038
41046-01	5	65	5	0.038
41047-01	5	65	10	0.038
41048-01	5	65	20	0.038
41049-01	5	100	5	0.038
41050-01	5	100	10	0.038
41051-01	5	100	20	0.038
41052-01	5	100	30	0.038
41053-01	5	100	40	0.038
41054-01	5	100	50	0.038
41055-01	5	135	5	0.038
41056-01	5	135	10	0.038
41057-01	5	135	20	0.038
41058-01	5	135	30	0.038
41059-01	5	135	40	0.038
41060-01	5	135	50	0.038

Snare



Amplatz Goose Neck™
Snare

Amplatz Goose Neck™

Snare Kit

The Amplatz Goose Neck snare kit is intended for use in the retrieval of atraumatic foreign bodies. Each kit contains one snare, one catheter, one introducer, and one torque device. The snare is constructed of nitinol cable and a gold-plated tungsten loop. Because of the snare's preformed loop and superelastic construction, it can be introduced through catheters without the risk of snare deformation.[†] The snare catheter has a platinum-iridium radiopaque marker band. The Amplatz Goose Neck snare kit is intended for the cardiovascular and peripheral vascular systems.

[†]Amplatz snare kit IFU.

Reference number	Loop diameter (mm)	Snare length (cm)	Catheter size (F)	Catheter length (cm)
GN500	5	120	4	102
GN1000	10	120	4	102
GN1001	10	65	4	48
GN1500	15	120	6	102
GN2000	20	120	6	102
GN2500	25	120	6	102
GN2501	25	65	6	48
GN3000	30	120	6	102
GN3500	35	120	6	102

Amplatz Goose Neck™

Microsnare Kit

The Amplatz Goose Neck microsnare kit is intended for the coronary and peripheral vascular systems and the extracranial neurovascular anatomy. Each kit contains one microsnare, one microsnare catheter, one introducer, and one torque device.

Reference number	Loop diameter (mm)	Snare length (cm)	Catheter size (F)	Catheter length (cm)
SK200	2	175	2.3-3	150
SK201	2	200	2.3-3	175
SK400	4	175	2.3-3	150
SK401	4	200	2.3-3	175
SK700	7	175	2.3-3	150
SK701	7	200	2.3-3	175

Guidewires



Wholey™*
Guidewire System



Nitrex™
Guidewire

Wholey™*

Guidewire System

The Wholey™* guidewire is intended to facilitate the placement and exchange of interventional devices during diagnostic or therapeutic interventional procedures. Wholey™* guidewires provide torqueability and lubricity, allowing interventionalists to approach challenging cases with confidence. Each Wholey™* wire, excluding the extension system, is packaged with a torque device.

Reference number	Description	Stiffness profile	Tip style	Size outer dia. (in)	Size length (cm)	Quantity
WWFS35145	Floppy tip, extension compatible	Floppy	Straight/Shapeable	0.035	145	3/Pkg
WWFS35175	Floppy tip, extension compatible	Floppy	Straight/Shapeable	0.035	175	3/Pkg
WWFS35260	Floppy tip, exchange length	Floppy	Straight/Shapeable	0.035	260	3/Pkg
WWFS35300	Floppy tip, exchange length	Floppy	Straight/Shapeable	0.035	300	3/Pkg
WWIJ35145	Modified J tip, extension compatible	Intermediate	Modified J/Shapeable	0.035	145	3/Pkg
WWIJ35175	Modified J tip, extension compatible	Intermediate	Modified J/Shapeable	0.035	175	3/Pkg
WWIJ35260	Modified J tip, exchange length	Intermediate	Modified J/Shapeable	0.035	260	3/Pkg
WWIJ35300	Modified J tip, exchange length	Intermediate	Modified J/Shapeable	0.035	300	3/Pkg
WWSS35145	Standard tip, extension compatible	Standard	Straight/Shapeable	0.035	145	3/Pkg
WWSS35175	Standard tip, extension compatible	Standard	Straight/Shapeable	0.035	175	3/Pkg
WWSS35260	Standard tip, exchange length	Standard	Straight/Shapeable	0.035	260	3/Pkg
WWSS35300	Standard tip, exchange length	Standard	Straight/Shapeable	0.035	300	3/Pkg
WWES35001	Extension system	Standard	Straight/Shapeable	0.035	155	3/Pkg
WWTD35001	Kendall™* torque device	NA	NA	0.025"-0.038"	NA	10/Pkg

Babywire™

Double-ended Nitinol Guidewire

The Babywire guidewire assists in the placement of initial catheters and in the exchange of catheters in small vessel anatomy. Babywire guidewires are straight 0.012" nitinol guidewires designed with double-ended round tips and flexible ends.

Reference number (10 package)	Diameter (in)	Length (cm)
BW1200	0.012	18
BW1201	0.012	50

Nitrex™

Guidewire

Nitrex guidewires are constructed of a superelastic nitinol core wire with a gold-plated tungsten coil for visibility. The wires have proprietary silicone coating for ease in placement. A torque device is packaged with 0.014" and 0.018" Nitrex guidewires. The 0.014" and 0.018" Nitrex guidewires are indicated for use in the peripheral and coronary vasculatures. The 0.025" and 0.035" Nitrex guidewires are indicated for use in the peripheral vasculature.

Reference number (3/Package) [†]	Diameter (in)	Length (cm)	Tip style	Tip length (cm)	Tip shape	Tip angle (°)
0.014"						
N140801	0.014	80	INT	5	A	15
N141802	0.014	180	INT	5	A	15
N143001	0.014	300	INT	5	A	15
0.018"						
N180601	0.018	60	INT	5	S	0
N180603	0.018	60	INT	7	S	0
N180801	0.018	80	STD	2	S	0
N180802	0.018	80	INT	5	A	15
N181804	0.018	180	STD	2	S	0
N181805	0.018	180	INT	5	A	15
N181806	0.018	180	FLOP	20	A	15
N183001	0.018	300	STD	2	S	0
N183002	0.018	300	INT	5	A	15
0.025"						
N251801	0.025	180	INT	8	A	15
N251802	0.025	180	STD	2	S	0
N252601	0.025	260	INT	8	A	15
0.035" Flexible Shaft						
N351451	0.035	145	INT	15	S	0
N351452	0.035	145	INT	15	A	45
N351803	0.035	180	INT	15	S	0
N352601	0.035	260	INT	15	A	45
N354001	0.035	400	INT	15	S	0
0.035" Stiff Shaft						
N350801	0.035	80	INT	9	S	0
N351453	0.035	145	FLOP	14	A	45
N351454	0.035	145	INT	9	S	0
N351455	0.035	145	FLOP	14	S	0
N351804	0.035	180	INT	9	S	0
N351805	0.035	180	STD	4	A	45
N352602	0.035	260	FLOP	14	S	0
N352603	0.035	260	STD	4	A	45
N352604	0.035	260	INT	9	S	0
N353001	0.035	300	INT	9	S	0
N354002	0.035	400	INT	9	S	0

[†]Torque devices included on 0.014" and 0.018" wire sizes.

Support Catheters and Guide Catheters



TrailBlazer™
Angled Support Catheter

TrailBlazer™

Support Catheter

The TrailBlazer catheter is an over-the-wire, single-lumen, seamless catheter with three embedded radiopaque markers, an atraumatic tapered tip, and a 40 cm hydrophilic distal tip coating. The TrailBlazer catheter is designed for high visibility, optimal wire support, and ease of lesion entry for difficult-to-cross lesions.

Reference number (5/package)	Guidewire compatibility	Working length (cm)	Minimum guide sheath (F)	Minimum introducer sheath (F)	Marker band spacing (mm)
SC-014-135	0.014	135	5	4	15
SC-014-150	0.014	150	5	4	15
SC-018-090	0.018	90	5	4	15
SC-018-135	0.018	135	5	4	15
SC-018-150	0.018	150	5	4	15
SC-035-065	0.035	65	6	5	50
SC-035-090	0.035	90	6	5	50
SC-035-135	0.035	135	6	5	50
SC-035-150	0.035	150	6	5	50

TrailBlazer™

Angled Support Catheter

With an angled tip and braided shaft design, the TrailBlazer angled support catheter has exceptional pushability and directionality to reach and cross lesions.

Reference number (5/Package)	Guidewire compatibility	Working length (cm)	Minimum guide sheath (F)	Minimum introducer sheath (F)	Marker and spacing (mm)
ASC-014-090	0.014	90	5	4	15
ASC-014-135	0.014	135	5	4	15
ASC-014-150	0.014	150	5	4	15
ASC-018-090	0.018	90	5	4	15
ASC-018-135	0.018	135	5	4	15
ASC-018-150	0.018	150	5	4	15
ASC-035-065	0.035	65	5	4	50
ASC-035-090	0.035	90	5	4	50
ASC-035-135	0.035	135	5	4	50
ASC-035-150	0.035	150	5	4	50

Launcher™

Peripheral Guide Catheter

The Launcher peripheral guide catheter is designed for multiple interventional approaches.



Reference number	French size (F)	Length (cm)	Curve style
Renal Curve			
LA6PK1W	6	47	PK1
LA7PK1W	7	47	PK1
LA8PK1W	8	47	PK1
Hockey Stick			
LA6MPHK	6	55	MPH
LA7MPHK	7	55	MPH
LA8MPHK	8	55	MPH
Renal Double Curve			
LA6RDCK	6	55	RDC
LA7RDCK	7	55	RDC
LA8RDCK	8	55	RDC
Sheperd's Crook			
LA6SCR40K	6	55	SCR
LA7SCR40K	7	55	SCR
LA8SCR40K	8	55	SCR

Reference number	French size (F)	Length (cm)	Curve style
Multipurpose			
LA6MP1K	6	55	MP1
LA7MP1K	7	55	MP1
LA8MP1K	8	55	MP1
Champ			
LA6CHAMP15K	6	55	Champ 1.5
LA7CHAMP15K	7	55	Champ 1.5
LA8CHAMP15K	8	55	Champ 1.5
LA6CHAMP20K	6	55	Champ 2.0
LA7CHAMP20K	7	55	Champ 2.0
LA8CHAMP20K	8	55	Champ 2.0
LA6CHAMP25K	6	55	Champ 2.5
LA7CHAMP25K	7	55	Champ 2.5
LA8CHAMP25K	8	55	Champ 2.5

Brief statements and indications for use



Abre™ venous self-expanding stent system

Intended Use/Indications: The Abre™ venous self-expanding stent system (Abre™ stent system) is indicated for use in the iliofemoral veins for the treatment of symptomatic venous outflow obstruction.

Contraindications: Do not use the Abre™ stent system with patients with known hypersensitivity to nickel titanium (nitinol), with patients who are judged to have a lesion that prevents complete inflation of a balloon dilatation catheter or proper placement of the stent or the stent delivery system, and with patients in whom anticoagulant or antiplatelet therapy is contraindicated.

Potential Adverse Effects of the Device on Health: The potential adverse effects (e.g., complications) associated with the use of the Abre™ stent system include, but are not limited to, access failure, access site infection, allergic reaction to contrast medium or procedure medications; aneurysm; AV fistula; bleeding; bruising; death; device breakage; device maldeployment; edema; embolization; fever; hematoma; hypertension; hypotension, nausea, or other vasovagal response; infection; myocardial infarction, arrhythmia, or other cardiovascular insufficiency; open surgical repair; pain; pseudoaneurysm; renal insufficiency or renal failure (new or worsening); respiratory distress or pulmonary embolism; sepsis; stent fracture; stent malapposition; stent malposition; stent migration; stroke, paradoxical embolism, transient ischemic attack, or intracerebral hemorrhage; tissue necrosis; venous occlusion, restenosis, or thrombosis, within or outside of stented segment; and vessel damage, including intimal injury, dissection, perforation, or rupture.

Admiral™ Xtreme™ 0.035" PTA Balloon Catheter

Indications for Use: The Admiral Xtreme PTA Balloon Dilatation Catheter is intended to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

Chocolate™ PTA Balloon

Indications for Use: The Chocolate™ PTA Balloon Catheter is intended for balloon dilatation of lesions in the peripheral vasculature, including the iliac, femoral, iliofemoral, popliteal, and infra-popliteal arteries.

Babywire™ Guidewire

Indications for Use: The Babywire™ Guidewire is intended for assisting in the placement of initial catheters or in the exchange of catheters in the small vessel anatomy. The Babywire guidewire is compatible with a 24 gauge needle or 0.7 mm (2 Fr) catheter.

ClosureFast™ Endovenous Radiofrequency (RFA) Ablation Catheter

Indications for Use: The ClosureFast™ endovenous radiofrequency ablation (RFA) catheter is intended for endovascular coagulation of blood vessels in patients with superficial vein reflux.

ClosureRFS™ Endovenous Radiofrequency Stylet

Indications for Use: The ClosureRFS™ stylet is intended for use in vessel and tissue coagulation, including treatment of incompetent (i.e., refluxing) perforator and tributary veins.

ClosureRFG™ Radiofrequency Generator

Indications for Use: The ClosureRFG generator is used with radiofrequency catheters intended for vessel and tissue coagulation.

Cragg-McNamara™ Valved Infusion Catheters

Indications for Use: The Cragg-McNamara™ Valved Infusion Catheter is intended to be used for the controlled selective infusion of physician-specified pharmacologic agents or radiopaque contrast media into the general vasculature. All pharmacologic agents utilized with the Micro Therapeutics Infusion Catheter should be fully prepared and used according to the instructions for use of the specific pharmacologic agent. The Micro Therapeutics Infusion Catheter is not intended for coronary, pediatric, or neonatal use.

Concerto™ detachable coil system

Indications for Use: The Concerto detachable coil system is indicated for arterial and venous embolizations in the peripheral vasculature.

Ellipsys™ System

The Ellipsys™ system is indicated for the creation of a proximal radial artery to perforating vein anastomosis via a retrograde venous access approach in patients with a minimum vessel diameter of 2.0 mm and less than 1.5 mm of separation between the artery and vein at the fistula creation site who have chronic kidney disease requiring dialysis.

Contraindications: The Ellipsys™ system is contraindicated for use in patients with target vessels that are < 2 mm in diameter. The Ellipsys™ System is contraindicated for use in patients who have a distance between the target artery and vein > 1.5 mm.

Warnings

- The Ellipsys™ system has only been studied for the creation of an AV fistula using the proximal radial artery and the adjacent perforating vein. It has not been studied in subjects who are candidates for surgical fistula creation at other locations, including sites distal to this location.
- The Ellipsys™ system is not intended to treat patients with significant vascular disease or calcification in the target vessels.
- The Ellipsys™ system has only been studied in subjects who had a patent palmar arch and no evidence of ulnar artery insufficiency.
- Use only with the Ellipsys™ Power Controller, Model No. AML-1001.
- The Ellipsys™ Catheter has been designed to be used with the 6 F Terumo Glidesheath Slender™. If using a different sheath, verify the catheter can be advanced through the sheath without resistance prior to use.

- Use ultrasound imaging to ensure proper placement of the catheter tip in the artery before retracting the sheath, since once the distal tip of the catheter has been advanced into the artery, it cannot be easily removed without creation of the anastomosis. If the distal tip is advanced into the artery at an improper location, complete the procedure and remove the catheter as indicated in the directions for use. It is recommended that a follow-up evaluation of the patient is performed using appropriate clinical standards of care for surgical fistulae to determine if any clinically significant flow develops that require further clinical action.

Precautions

- This product is sterilized by ethylene oxide gas.
 - Additional procedures are expected to be required to increase and direct blood flow into the AVF target outflow vein and to maintain patency of the AVF. Care should be taken to proactively plan for any fistula maturation procedures when using the device.
 - In the Ellipsys™ study, 99% of subjects required balloon dilatation (PTA) to increase flow to the optimal access vessel and 62% of subjects required embolization coil placement in competing veins to direct blood flow to the optimal access vessel. Prior to the procedure, care should be taken to assess the optimal access vessel for maturation, the additional procedures that may be required to successfully achieve maturation, and appropriate patient follow-up. Please refer to the "Arteriovenous Fistula (AVF) Maturation" section of the labeling for guidance about fistula flow, embolization coil placement, and other procedures to assist fistula maturation and maintenance.
 - The Ellipsys™ System is intended to only be used by physicians trained in ultrasound guided percutaneous endovascular interventional techniques using appropriate clinical standards for care for fistula maintenance and maturation including balloon dilatation and coil embolization.
 - Precautions to prevent or reduce acute or longer-term clotting potential should be considered. Physician experience and discretion will determine the appropriate anticoagulant/antiplatelet therapy for each patient using appropriate clinical standards of care.
- Potential Adverse Events:** Potential complications that may be associated with creation and maintenance of an arteriovenous fistula include, but may not be limited to, the following:
- Total occlusion, partial occlusion or stenosis of the anastomosis or

Brief statements and indications for use



- adjacent outflow vein
- Stenosis of the central AVF outflow requiring treatment per the treatment center's standard of care
- Failure to achieve fistula maturation
- Incomplete vessel ligation when using embolization coil to direct flow
- Steal Syndrome
- Hematoma
- Infection or other complications
- Need for vessel superficialization or other maturation assistance procedures.

Enteer™ Re-entry Catheter

Indications for Use: The Enteer™ Re-entry Catheter is indicated for directing, steering, controlling, and supporting a guidewire in order to access discrete regions of the peripheral vasculature. When used as part of the Peripheral System, the Enteer Catheter is indicated for use to facilitate the intraluminal placement of conventional guidewires beyond stenotic peripheral lesions (including chronic total occlusions) prior to placement of other interventional devices.

Enteer™ Re-entry Guidewire

Indications for Use: The Enteer™ Re-entry Guidewire is intended to facilitate placement of balloon dilatation catheters or other intravascular devices during percutaneous transluminal angioplasty (PTA). The Enteer Guidewire is not to be used in cerebral blood vessels. When used as part of the Peripheral System, the Enteer Guidewire is indicated for use to facilitate the intraluminal placement of conventional guidewires beyond stenotic peripheral lesions (including chronic total occlusions) prior to placement of other interventional devices.

EverCross™ 0.035" PTA Balloon Catheter

Indications for Use: The EverCross™ 0.035" OTW PTA Dilatation Catheter is intended to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, infrapopliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

EverFlex™ Self-expanding Peripheral Stent System

Indication: The EverFlex™ Self-Expanding Peripheral Stent System is intended to improve luminal diameter in the treatment of symptomatic de novo or restenotic lesions up to 180 mm in length in the native superficial femoral artery and/or proximal popliteal arteries with reference vessel diameters ranging from 4.5 mm-7.5 mm. The EverFlex Self-Expanding Peripheral Stent System is indicated for improving luminal diameter in patients with atherosclerotic disease of the common and/or external iliac arteries up to and including 100 mm in length, with a reference vessel diameter of 4.5 mm-7.5 mm. The Protégé EverFlex Self-expanding Biliary Stent System is intended as a palliative treatment of malignant neoplasms in the biliary tree.

Contraindications: Use of the EverFlex™ Self-Expanding Peripheral Stent System is contraindicated in patients with known hypersensitivity to nickel titanium and in patients contraindicated for anticoagulant and/or antiplatelet therapy, patients who are judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or stent delivery system.

Potential Adverse Events: Potential adverse events which may be associated with the use of a stent in the SFA, proximal popliteal arteries and iliac arteries include, but are not limited to: abrupt or sub-acute closure, allergic reaction to device materials or procedure medications, allergic reaction to nitinol, amputation, aneurysm, angina, arrhythmia, arterio-venous fistula, artery injury (e.g., dissection, perforation, or rupture), bleeding requiring transfusion, bruising, contrast medium reaction/renal failure, death, device breakage, edema, embolism, failure to deploy stent, fever, gastrointestinal bleeding due to anticoagulation, hematoma, hypertension/hypotension, infection, inflammation, intraluminal thrombus, myocardial infarction, pain, partial

stent deployment, pseudoaneurysm, renal failure, renal insufficiency, restenosis, sepsis, shock, stent collapse or fracture, stent migration, stent misplacement, stroke, surgical or endovascular intervention, thrombosis/occlusion of the stent transient ischemic attack, venous thromboembolism, vessel spasm, and worsening claudication or rest pain. Potential adverse events which may be associated with the use of a stent in the biliary include, but are not limited to: infection secondary to contamination of the stent may lead to cholangitis, hemobilia, peritonitis, or abscess; the stent may migrate from the site of implant down the biliary tract, overstretching the duct may result in rupture, persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this implant, device breakage, failure to deploy the stent, partial stent deployment, stent collapse or fracture, stent misplacement, and surgical intervention.

EverFlex™ Self-expanding Peripheral Stent with Entrust™ Delivery System

Indication: The EverFlex™ Self-Expanding Peripheral Stent with Entrust™ Delivery System is intended to improve luminal diameter in the treatment of symptomatic de-novo or restenotic lesions up to 140 mm in length in the native Superficial Femoral Artery (SFA) and/or proximal popliteal arteries with reference vessel diameters ranging from 4.5 mm-7.5 mm.

Contraindications: Use of the EverFlex™ Self-Expanding Peripheral Stent with Entrust™ Delivery System is contraindicated in patients with known hypersensitivity to nickel titanium; patients contraindicated for anticoagulant and/or antiplatelet therapy; patients who have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or stent delivery system. The EverFlex™ Self-Expanding Peripheral Stent with Entrust™ Delivery System is contraindicated for use in the carotid artery.

Potential adverse events which may be associated with the use of a stent in the SFA and proximal popliteal arteries include, but are not limited to: abrupt or subacute closure, allergic reaction to device materials or procedure medications, allergic reaction to nitinol, amputation, aneurysm, angina, arrhythmia, arteriovenous fistula, artery perforation or rupture, bleeding requiring transfusion, bruising, contrast medium reaction /renal failure, death, device breakage, dissection or intimal flap, edema, embolism, failure to deploy stent, fever, gastrointestinal bleeding due to anticoagulation, hematoma, hypertension/hypotension, infection, inflammation, intraluminal thrombus, myocardial infarction, pain, partial stent deployment, pseudoaneurysm, renal failure requiring dialysis, renal insufficiency (new or worsening), restenosis, sepsis, shock, stent collapse or fracture, stent migration, stent misplacement, stroke, surgical or endovascular intervention, thrombosis/occlusion of the stent, transient ischemic attack, venous thromboembolism, and vessel spasm.

Fortrex™ HP PTA Balloon Catheter

Indications for Use: The Fortrex™ HP PTA balloon catheter is intended to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

Goose Neck™ Microsnare

Indications for Use: The Amplatz Goose Neck microsnare kit is intended for use in the retrieval and manipulation of atraumatic foreign bodies located in the coronary and peripheral cardiovascular system and the extracranial neurovascular anatomy.

Goose Neck™ Snare

Indications for Use: The Amplatz Goose Neck snare is intended for use in the cardiovascular system or hollow viscus to retrieve and manipulate foreign objects. Manipulation procedures include indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping, and central venous access venipuncture procedure assistance.

Brief statements and indications for use



HawkOne™ Directional Atherectomy System

Indications for Use: The HawkOne directional atherectomy system is intended for use in atherectomy of the peripheral vasculature. The HawkOne catheter is indicated for use in conjunction with the SpiderFX™ embolic protection device in the treatment of severely calcified lesions. The HawkOne catheter is NOT intended for use in the coronary, carotid, iliac or renal vasculature.

IN.PACT™ 018 Drug-Coated Balloon and IN.PACT™ Admiral™ Drug-Coated Balloon

Indications for Use: The IN.PACT Admiral Paclitaxel-coated PTA Balloon Catheter and IN.PACT 018 Paclitaxel-coated PTA Balloon Catheter are indicated for percutaneous transluminal angioplasty, after appropriate vessel preparation, of de novo, restenotic, or in-stent restenotic lesions with lengths up to 360 mm in superficial femoral or popliteal arteries with reference vessel diameters of 4-7 mm.

Contraindications:

The IN.PACT Admiral DCB and IN.PACT 018 DCB are contraindicated for use in: Coronary arteries, renal arteries, and supra-aortic/cerebrovascular arteries,

Patients who cannot receive recommended antiplatelet and/or anticoagulant therapy

Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the delivery system

Patients with known allergies or sensitivities to paclitaxel

Women who are breastfeeding, pregnant, or are intending to become pregnant or men intending to father children. It is unknown whether paclitaxel will be excreted in human milk and whether there is a potential for adverse reaction in nursing infants from paclitaxel exposure.

Warnings:

• **A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients.**

• Use the product prior to the Use-by Date specified on the package. Contents are supplied sterile. Do not use the product if the inner packaging is damaged or opened.

• Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium (equal parts contrast medium and saline solution).

• Do not exceed the rated burst pressure (RBP). The RBP is based on the results of in vitro testing. Use of pressures higher than RBP may result in a ruptured balloon with possible intimal damage and dissection. [IN.PACT Admiral DCB: The RBP is 14 atm (1419 kPa) for all balloons except the 200 and 250 mm balloons. For the 200 and 250 mm balloons the RBP is 11 atm (1115 kPa) - IN.PACT 018 DCB: The RBP is 10 atm (1013 kPa) for all balloons.]

• The safety and effectiveness of using multiple IN.PACT Admiral DCB, or multiple IN.PACT 018 DCB, with a total drug dosage exceeding 34,854 µg of paclitaxel in a patient has not been clinically evaluated.

Precautions

The safety and effectiveness of the IN.PACT Admiral DCB (0.035 in guidewire compatible), as established in the clinical studies that were performed primarily via femoral access, can be considered supportive for the IN.PACT 018 DCB. Vessel preparation using only pre-dilatation was studied in the IN.PACT Admiral DCB clinical studies. Other methods of vessel preparation, such as atherectomy, have not been studied clinically. The IN.PACT 018 DCB has not been evaluated in a clinical study.

The IN.PACT Admiral DCB and IN.PACT 018 DCB should only be used

by physicians trained in percutaneous transluminal angioplasty (PTA).

The IN.PACT Admiral DCB and IN.PACT 018 DCB are designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Assess risks and benefits before treating patients with a history of severe reaction to contrast agents.

The safety and effectiveness of the IN.PACT Admiral DCB or IN.PACT 018 DCB used in conjunction with other drug-eluting stents or drug-coated balloons in the same procedure or following treatment failure has not been evaluated.

The extent of the patient's exposure to the drug coating is directly related to the number of balloons used. Refer to the Instructions for Use (IFU) for details regarding the use of multiple balloons and paclitaxel content.

The use of the IN.PACT Admiral DCB and IN.PACT 018 DCB carries the risks associated with percutaneous transluminal angioplasty, including thrombosis, vascular complications, and/or bleeding events.

The IN.PACT Admiral DCB and IN.PACT 018 DCB are not intended for the expansion or delivery of a stent.

Potential Adverse Effects

The potential adverse effects (e.g., complications) associated with the use of the device are: abrupt vessel closure; access site pain; allergic reaction to contrast medium, antiplatelet therapy, or catheter system components (materials, drugs, and excipients); amputation/loss of limb; arrhythmias; arterial aneurysm; arterial thrombosis; arteriovenous (AV) fistula; death; dissection; embolization; fever; hematoma; hemorrhage; hypotension/hypertension; inflammation; ischemia or infarction of tissue/ organ; local infection at access site; local or distal embolic events; perforation or rupture of the artery; pseudoaneurysm; renal insufficiency or failure; restenosis of the dilated artery; sepsis or systemic infection; shock; stroke; systemic embolization; vessel spasms or recoil; vessel trauma which requires surgical repair.

Potential complications of peripheral balloon catheterization include, but are not limited to: balloon rupture; detachment of a component of the balloon and/or catheter system; failure of the balloon to perform as intended; failure to cross the lesion. Although systemic effects are not anticipated, potential adverse events that may be unique to the paclitaxel drug coating include, but are not limited to: allergic/immunologic reaction; alopecia; anemia; gastrointestinal symptoms; hematologic dyscrasia (including leucopenia, neutropenia, thrombocytopenia); hepatic enzyme changes; histologic changes in vessel wall, including inflammation, cellular damage, or necrosis; myalgia/arthralgia; myelosuppression; peripheral neuropathy. Refer to the Physician's Desk Reference for more information on the potential adverse effects observed with paclitaxel. There may be other potential adverse effects that are unforeseen at this time. Please reference appropriate product Instructions for Use for a detailed list of indications, warnings, precautions, and potential adverse effects. This content is available electronically at manuals.medtronic.com.

IntraStent™ Biliary Stent System

Indications for Use: The stent is intended as a palliative treatment of malignant neoplasms in the biliary tree.

WARNING: The safety and effectiveness of this device for use in the vascular system have not been established.

Launcher™ Peripheral Guide Catheter

Indications for Use: The Launcher™ Peripheral Guide Catheter is designed to provide a pathway through which therapeutic devices are introduced. The guiding catheter is intended to be used in the coronary or peripheral vascular system.

IN.PACT™ AV Paclitaxel-coated PTA Balloon

INDICATIONS FOR USE: The IN.PACT™ AV Paclitaxel-coated

Brief statements and indications for use



PTA Balloon Catheter is indicated for percutaneous transluminal angioplasty, after appropriate vessel preparation, for the treatment of obstructive lesions up to 100 mm in length in the native arteriovenous dialysis fistulae with reference vessel diameters of 4 to 12 mm.

CONTRAINDICATIONS:

The IN.PACT AV DCB is contraindicated for use in the following anatomy and patient types:

- Coronary arteries, renal arteries, and supra-aortic/cerebrovascular arteries
- Patients who cannot receive recommended antiplatelet and/or anticoagulant therapy
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the delivery system
- Patients with known allergies or sensitivities to paclitaxel
- Women who are breastfeeding, pregnant, or are intending to become pregnant, or men intending to father children. It is unknown whether paclitaxel will be excreted in human milk and whether there is a potential for adverse reaction in nursing infants from paclitaxel exposure

WARNINGS:

A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Inadequate information is available to evaluate the potential mortality risk associated with the use of paclitaxel-coated devices for the treatment of other diseases/conditions, including this device indicated for use in arteriovenous dialysis fistulae. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options for their specific disease/condition with their patients.

- Use the product prior to the Use-by date specified on the package.
- Contents are supplied sterile. Do not use the product if the inner packaging is damaged or opened.
- Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium (equal parts contrast medium and saline solution).
- Do not move the guidewire during inflation of the IN.PACT AV DCB.
- Do not exceed the rated burst pressure (RBP). The RBP is based on the results of in vitro testing. Use of pressures higher than RBP may result in a ruptured balloon with possible intimal damage and dissection.
- The safety of using multiple IN.PACT AV DCBs with a total drug dosage exceeding 15,105 µg paclitaxel has not been evaluated clinically.

PRECAUTIONS:

- This product should only be used by physicians trained in percutaneous transluminal angioplasty (PTA).
- Assess risks and benefits before treating patients with a history of severe reaction to contrast agents. Identify allergic reactions to contrast media and antiplatelet therapy before treatment and consider alternatives for appropriate management prior to the procedure. This product is not intended for the expansion or delivery of a stent.
- Do not use the IN.PACT AV DCB for pre-dilatation or for post-dilatation.
- This product is designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- The use of this product carries the risks associated with percutaneous transluminal angioplasty, including thrombosis, vascular complications, and/or bleeding events
- The safety and effectiveness of the IN.PACT AV DCB used in conjunction with other drug-eluting stents or drug-coated balloons in the same procedure has not been evaluated.
- The extent of the patient's exposure to the drug coating is directly related to the number of balloons used. Refer to the Instructions for Use (IFU) for details regarding the use of multiple balloons and

paclitaxel content.

- Appropriate vessel preparation, as determined by the physician to achieve residual stenosis of $\leq 30\%$, is required prior to use of the IN.PACT AV DCB. Vessel preparation of the target lesion using high-pressure PTA for pre-dilatation was studied in the IN.PACT AV Access clinical study. Other methods of vessel preparation, such as atherectomy, have not been studied clinically with IN.PACT AV DCB.

POTENTIAL ADVERSE EFFECTS: Potential adverse effects which may be associated with balloon catheterization may include, but are not limited to, the following: abrupt vessel closure, allergic reaction, arrhythmias, arterial or venous aneurysm, arterial or venous thrombosis, death, dissection, embolization, hematoma, hemorrhage, hypotension/hypertension, infection, ischemia or infarction of tissue/organ, loss of permanent access, pain, perforation or rupture of the artery or vein, pseudoaneurysm, restenosis of the dilated vessel, shock, stroke, vessel spasms or recoil.

Potential complications of peripheral balloon catheterization include, but are not limited to, the following: balloon rupture, detachment of a component of the balloon and/or catheter system, failure of the balloon to perform as intended, failure to cross the lesion. These complications may result in adverse effects.

Although systemic effects are not anticipated, potential adverse effects not captured above that may be unique to the paclitaxel drug coating include, but are not limited to, the following: allergic/immunologic reaction, alopecia, anemia, gastrointestinal symptoms, hematologic dyscrasia (including leucopenia, neutropenia, thrombocytopenia), hepatic enzyme changes, histologic changes in vessel wall, including inflammation, cellular damage, or necrosis, myalgia/arthralgia, myelosuppression, peripheral neuropathy.

Refer to the Physician's Desk Reference for more information on the potential adverse effects observed with paclitaxel. There may be other potential adverse effects that are unforeseen at this time.

Please reference appropriate product Instructions for Use for a detailed list of indications, warnings, precautions and potential adverse effects. This content is available electronically at www.manuals.medtronic.com.

Mo.Ma™ Ultra Cerebral Protection Device

Indications for Use: The Mo.Ma™ Ultra Proximal Cerebral Protection Device is indicated as an embolic protection system to contain and remove embolic material (thrombus/debris) while performing angioplasty and stenting procedures involving lesions of the internal carotid artery and/or the carotid bifurcation.

The reference diameter of the external carotid artery should be between 3-6 mm and the reference diameter of the common carotid artery should be between 5-13 mm.

MVP™ micro vascular plug system

Indications for Use: The MVP micro vascular plug system is indicated to obstruct or reduce the rate of blood flow in the peripheral vasculature.

NanoCross™ Elite 0.014" PTA Balloon Catheter

Indications for Use: The NanoCross™ Elite 0.014" Over-the-Wire PTA Balloon Dilatation Catheter is intended to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

Nitrex™ Guidewire

Indications for Use: The 0.014 in. (0.36 mm) and 0.018 in. (0.46 mm) diameter NITREX Nitinol Guidewires are intended for use in the peripheral and coronary vasculature. The 0.025 in. (0.64 mm) and 0.035 in. (0.89 mm) diameter NITREX Nitinol Guidewires are indicated for use in the peripheral vasculature.

Pacific™ Plus 0.018" PTA Balloon Catheter

Brief statements and indications for use



Indications for Use: The Pacific™ Plus PTA catheter is intended to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, infrapopliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. Catheters with balloon diameters of 4-7 mm (in all lengths) are also indicated for stent postdilatation in the peripheral arteries.

Protégé™ GPS™ Self-expanding Peripheral and Biliary Stent System

Indication: The Protégé™ GPS™ Self-Expanding Peripheral Stent Systems is indicated for improving luminal diameter in patients with atherosclerotic disease of the common and/or external iliac arteries up to and including 100 mm in length, with reference vessel diameters of 7.5-11 mm. The stent is intended as a palliative treatment of malignant neoplasms in the biliary tree.

Contraindications: Use of the Protégé™ GPS™ Self-Expanding Peripheral Stent System is contraindicated in patients with known hypersensitivity to nickel titanium; patients contraindicated for anticoagulant and/or antiplatelet therapy; patients who have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or stent delivery system.

Potential Adverse Events: Potential adverse events which may be associated with the use of a stent in the common and/or external iliac arteries include, but are not limited to: Abrupt or sub-acute closure, Allergic reaction to device materials or procedure medications, Allergic reaction to Nitinol, Amputation, Aneurysm, Angina, Arrhythmia, Arterio-venous fistula, Artery injury (e.g., dissection, perforation, or rupture), Bleeding requiring transfusion, Bruising, Contrast medium reaction/renal Failure, Death, Device breakage, Edema, Embolism, Failure to deploy stent, Fever, Gastrointestinal bleeding due to Anticoagulation, Hematoma, Hypertension/Hypotension, Infection, Inflammation, Intraluminal thrombus, Myocardial infarction Pain, Partial stent deployment, Pseudoaneurysm, Renal failure, Renal insufficiency, Restenosis, Sepsis, Shock, Stent collapse or fracture, Stent migration, Stent misplacement, Stroke, Surgical or endovascular Intervention, Thrombosis/occlusion of the stent, Transient ischemic attack, Venous thromboembolism, Vessel spasm, Worsening claudication or rest, pain. Potential adverse events which may be associated with the use of a stent in the biliary include, but are not limited to: infection secondary to contamination of the stent may lead to cholangitis, hemobilia, peritonitis, or abscess; the stent may migrate from the site of implant down the biliary tract, overstretching the duct may result in rupture, persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this implant, device breakage, failure to deploy the stent, partial stent deployment, stent collapse or fracture, stent misplacement, and surgical intervention.

Protégé™ RX Self-expanding Carotid Stent System

Indications: The Protégé™ RX Carotid Stent System, when used in conjunction with the ev3 embolic protection system, is indicated for the treatment of patients at high risk for adverse events from carotid endarterectomy who require percutaneous carotid revascularization and meet the following criteria: 1. Patients with carotid artery stenosis ($\geq 50\%$ for symptomatic patients by ultrasound or angiography or $\geq 80\%$ for asymptomatic patients by ultrasound or angiography) of the Common or Internal Carotid Artery, AND 2. Patients must have a reference vessel diameter within the range of 4.5 mm and 9.5 mm at the target lesion.

Contraindications: Use of the Protégé RX Carotid Stent System is contraindicated under these circumstances: Patients in whom anticoagulant, antiplatelet therapy or thrombolytic drugs is contraindicated; patients with vascular tortuosity or anatomy, which precludes the safe introduction of the sheath, guide catheter, embolic protection system, or stent system; patients with known hypersensitivity to nickel-titanium; patients with uncorrected bleeding disorders; lesions in the ostium of the common carotid artery.

WARNING: Only physicians who have received appropriate training and are familiar with the principles, clinical applications, complications, side effects and hazards commonly associated with carotid interventional

procedures should use this device.

- The device is provided sterile for single use only. DO NOT reprocess or resterilize. Reprocessing and resterilizing increase the risk of patient infection and risk of compromised device performance.
- Do not use the product after the "Use By Date" printed on the package.
- Maintain the patient's Activated Clotting Time (ACT) at > 250 seconds throughout PROTÉGÉ RX Carotid Stent System usage to prevent thrombus formation on the device.
- Maintain continuous flush while removing and reinserting devices on the guidewire. Perform all exchanges slowly to prevent air embolism or trauma to the artery.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stent and may cause acute closure of the vessel, requiring additional intervention (carotid endarterectomy, further dilatation, or placement of additional stents).
- The stent may cause thrombus, distal embolization or may migrate from the site of implant down the arterial lumen. Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration. In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infection, pseudoaneurysm or fistulization, surgical removal of the stent may be required.
- Overstretching of the artery may result in rupture and life-threatening bleeding.
- Never withdraw or move an intravascular device against any resistance until the cause is determined. Advancing with such resistance may lead to embolization of debris, and vessel and/or device damage.
- Frequently observe the PROTÉGÉ RX Carotid Stent System under fluoroscopy during stent deployment.
- Exercise caution when advancing or withdrawing the PROTÉGÉ RX Carotid Stent System through any previously placed devices.
- Allow for and maintain adequate distance between the embolic protection device and the stent delivery system or deployed stent to avoid potential entanglement.
- Ensure optimal positioning of the stent prior to deployment. Once deployment is initiated, the stent cannot be repositioned or recaptured. Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the carotid vasculature and/or the vascular access site. Complications may include death, stroke, bleeding, hematoma or pseudoaneurysm.

Potential Adverse Events: Potential adverse events which may be associated with the use of a stent in the carotid arteries include, but are not limited to: Abrupt closure, Allergic reactions to procedural medications, contrast dye or device materials, Amaurosis fugax, Aneurysm, Angina/coronary ischemia, Arrhythmia, Arterial occlusion or thrombosis at puncture site or remote site, Arteriovenous fistula, Bacteremia or septicemia, Bleeding from anticoagulant or antiplatelet medications, Bleeding, with or without transfusion, Cerebral edema, Cerebral hemorrhage, Cerebral ischemia or transient ischemic attack (TIA), Congestive heart failure (CHF), Death, Detachment of a component of the device system, Embolism (air, tissue, thrombus), Emergent or urgent endarterectomy surgery (CEA), Fever, Filter thrombosis or occlusion, Fluid overload, Groin hematoma, with or without surgical repair, Hemorrhage, with or without transfusion, Hyperperfusion syndrome, Hypotension or hypertension, Infection and/or pain at the puncture site, Ischemia or infarction of tissue/organ, Myocardial infarction (MI), Pain (head, neck), Pseudoaneurysm (femoral), Renal failure/insufficiency (new or worsening), Restenosis of stented segment, Seizure, Severe unilateral headache, Slow/no flow during procedure, Stent/filter collapse or fracture, Stent/filter entanglement or damage, Stent/filter failure to deploy, Stent embolization, migration, or misplacement, Stent or vessel thrombosis/occlusion, Stroke/cerebrovascular accident (CVA), Total occlusion of carotid artery, Vessel dissection, flap, perforation, or rupture, Vessel spasm or recoil.

RapidCross™ 0.014" Rapid Exchange PTA Balloon Catheter

Indications for Use: The RapidCross™ PTA Rapid Exchange Balloon

Brief statements and indications for use



Dilatation Catheter is intended to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

Rebar™ Micro Catheter

Indications for Use: The Rebar micro catheter is intended for the controlled selective infusion of physician-specified therapeutic agents or contrast media into the vasculature of the peripheral and neuro anatomy.

SilverHawk™ Plaque Excision System

Indications for Use: The SilverHawk peripheral plaque excision system is intended for use in atherectomy of the peripheral vasculature. The catheter is NOT intended for use in the coronary, carotid, iliac or renal vasculature.

SpiderFX™ Embolic Protection Device

Indications for Use: Lower Extremity (LE) Interventions

The SpiderFX™ Embolic Protection Device is indicated for use as a guidewire and embolic protection system to contain and remove embolic material in conjunction with the TurboHawk™ Peripheral Plaque Excision System, either during standalone procedures or together with PTA and/or stenting, in the treatment of severely calcified lesions in arteries of the lower extremities. The vessel diameter at the filter basket placement site should be between 3.0 mm and 6.0 mm.

Carotid Interventions

The SpiderFX Embolic Protection Device is indicated for use as a guidewire and embolic protection system to contain and remove embolic material (thrombus/debris) while performing angioplasty and stenting procedures in carotid arteries. The diameter of the artery at the site of filter basket placement should be between 3.0 mm and 7.0 mm.

Saphenous Vein Graft (SVG) Interventions

The SpiderFX Embolic Protection Device is indicated for use as an embolic protection system to contain and remove embolic material (thrombus/debris). The device also acts as the guidewire while performing percutaneous transluminal coronary angioplasty or stenting procedures in coronary saphenous vein bypass grafts with reference vessel diameters of 3.0 mm to 6.0 mm. The safety and effectiveness of this device as an embolic protection system has not been established in the cerebral vasculature.

TrailBlazer™ Support Catheter

Indications for Use: TrailBlazer™ Support Catheters are percutaneous, single-lumen catheters designed for use in the peripheral vascular system. TrailBlazer™ Support Catheters are intended to guide and support a guidewire during access of the vasculature, allow for wire exchanges, and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

TrailBlazer™ Angled Support Catheter

Indications for Use: TrailBlazer™ Angled Support Catheters are percutaneous, single-lumen catheters designed for use in the peripheral vascular system. TrailBlazer™ Support Catheters are intended to guide and support a guide-wire during access of the vasculature, allow for wire exchanges, and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

VenaSeal™ Closure System

Intended Use/Indications: The VenaSeal™ closure system (VenaSeal™ system) is indicated for use in the permanent closure of lower extremity superficial truncal veins, such as the great saphenous vein (GSV), through endovascular embolization with coaptation. The VenaSeal system is intended for use in adults with clinically symptomatic venous reflux as diagnosed by duplex ultrasound (DUS).

Contraindications: Separate use of the individual components of the VenaSeal closure system is contraindicated. These components must

be used as a system. The use of the VenaSeal system is contraindicated when any of the following conditions exist: previous hypersensitivity reactions to the VenaSeal™ adhesive or cyanoacrylates, acute superficial thrombophlebitis, thrombophlebitis migrans, acute sepsis. Potential Adverse Effects of the Device on Health: The potential adverse effects (e.g., complications) associated with the use of the VenaSeal system include, but are not limited to, adverse reactions to a foreign body (including, but not limited to, nonspecific mild inflammation of the cutaneous and subcutaneous tissue), arteriovenous fistula, bleeding from the access site, deep vein thrombosis (DVT), edema in the treated leg, embolization, including pulmonary embolism (PE), hematoma, hyperpigmentation, hypersensitivity or allergic reactions to cyanoacrylates, such as urticaria, shortness of breath, and anaphylactic shock, infection at the access site, pain, paresthesia, phlebitis, superficial thrombophlebitis, urticaria, erythema, or ulceration may occur at the injection site, vascular rupture and perforation, visible scarring. Warnings, precautions, and instructions for use can be found in the product labeling at <http://manuals.medtronic.com>.

Viance™ Crossing Catheter

Indications for Use: The Viance™ Catheter is intended for use with a guidewire to access discrete regions of the peripheral vasculature. When used as part of the Peripheral System, the Viance Catheter is indicated for use to facilitate the intraluminal placement of conventional guidewires beyond stenotic peripheral lesions (including chronic total occlusions) prior to placement of other interventional devices.

Visi-Pro™ Balloon-expandable Peripheral Stent System

Indications: The Visi-Pro™ Balloon-expandable Peripheral Stent System is indicated for improving luminal diameter in patients with atherosclerotic disease of the common and/or external iliac arteries up to 100 mm in length, with a reference vessel diameter of 5 to 10 mm. The Visi-Pro™ Balloon-expandable Biliary Stent System is intended as a palliative treatment of malignant neoplasms in the biliary tree.

Contraindications: Use of the Visi-Pro™ Balloon-expandable Peripheral Stent System is contraindicated in patients with known allergies to stainless steel or its components (for example nickel); patients contraindicated for anticoagulant and/or antiplatelet therapy; patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site; perforation at the angioplasty site evidenced by extravasation of contrast medium; aneurysm of the artery to be treated. All of the customary contraindications for PTA.

Potential Adverse Events: Potential adverse events which may be associated with the use of a stent in the iliac arteries include, but are not limited to: Abrupt or sub-acute closure, Allergic reaction to 316L stainless steel, Allergic reaction to device materials or procedure medications. Amputation, Aneurysm, Angina, Arrhythmia, Arterio-venous fistula, Artery injury (e.g., dissection, perforation or rupture), Bleeding requiring transfusion, Contrast medium reaction/renal failure, Death Device breakage, Embolism, Failure to deploy stent, Fever, Gastrointestinal bleeding due to anticoagulation, Hematoma, Hypertension/Hypotension, Infection, Inflammation, Intraluminal thrombus, Myocardial infarction, Pain, Partial stent deployment, Pseudoaneurysm, Renal insufficiency, Restenosis, Sepsis, Shock, Stent collapse or fracture, Stent migration, Stent misplacement, Stroke, Surgical or endovascular intervention, Thrombosis/occlusion of the stent, Transient increase in glomerular filtration rate, Transient ischemic attack, Venous thromboembolism, Vessel spasm, Worsening claudication or rest pain. Potential adverse events which may be associated with the use of this device in the biliary include, but are not limited to: infection secondary to contamination of the stent may lead to cholangitis, hemobilia, peritonitis, or abscess; the stent may migrate from the site of implant down the biliary tract, overstretching the duct may result in rupture, and because the stent is made of 316L stainless steel, persons with known allergic reactions to 316L stainless steel may suffer an allergic response to this implant.

Brief statements and indications for use



Wholey™ Guidewire System

Indications for Use: The Wholey™ guidewire system is intended to facilitate the placement and exchange of interventional devices during diagnostic or therapeutic interventional procedures. The guidewire can be torqued to facilitate navigation through tortuous arteries and/or avoid unwanted side branches.



Medtronic

[medtronic.com/PVH](https://www.medtronic.com/PVH)

UC202306641b EN ©2023 Medtronic. Medtronic and the Medtronic logo are trademarks of Medtronic.™*Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. Wholey is a trademark of Mark H. Wholey used under license. For distribution in the USA only. 07/2023